lumiraDx™

lumiraDx SARS-CoV-2 & RSV Quality Controls

For use with the LumiraDx Platform
For Professional Use Only
For *In Vitro* Diagnostic Use **IVD**

Product Name	REF
SARS-CoV-2 & RSV Quality Controls 2pk Carton - CE (EN/FR/DE/IT/ES/PT/NO/NL/PL/DA/PF/UK/PT/BR/PL/GR/RU/CZ/RU/UA)	L02108010002
SARS-CoV-2 & RSV Quality Controls 2pk Carton - CE (PT/BR/RSV/CZ/RU/UA)	L02108010002

SPEC-34881 R2	ART-01682 R2	Date of Rev 2022-05
---------------	--------------	---------------------

The LumiraDx SARS-CoV-2 & RSV Quality Controls (hereafter referred to as Quality Controls) are liquid quality control samples used with the LumiraDx Instrument (hereafter referred to as the Instrument) and the LumiraDx SARS-CoV-2 & RSV Test (hereafter referred to as SARS-CoV-2 & RSV Test). Read these instructions thoroughly before using the Quality Controls.

Inspect the Quality Controls packaging and contents for damage before use. Report any damage to LumiraDx Customer Services and do not use the kit if any damage is observed to the contents. To ensure that you are using the Instrument, the SARS-CoV-2 & RSV Test and the Quality Controls correctly, read the appropriate Platform User Manual, SARS-CoV-2 & RSV Test Product Insert and this entire product insert. In addition, please watch the LumiraDx Platform Training video available at lumiraDx.com.

Intended Use: The LumiraDx SARS-CoV-2 & RSV Quality Controls are intended for use by laboratory professionals/technicians/professionals for automated liquid quality control testing performed on the LumiraDx Instrument when used with the LumiraDx SARS-CoV-2 & RSV Test Strip. The Quality Controls provide users with assurance that the device is performing within specification. The LumiraDx SARS-CoV-2 & RSV Quality Control kit should be used whenever local testing guidelines require a QC to be run. The Quality Controls are intended for professional use only.

Summary and Explanation of the Test: The LumiraDx SARS-CoV-2 & RSV Quality Controls are liquid quality controls. The controls are specifically formulated and manufactured to ensure performance of the test and are used to verify the user's ability to properly perform the test and interpret the results.

It is the responsibility of each laboratory or healthcare setting using the LumiraDx SARS-CoV-2 & RSV Test to establish an adequate quality assurance program to ensure the performance of the device under its specific locations and conditions of use. Quality control requirements should be followed in conformance with local regulations or accreditation requirements and the user laboratory's standard quality control procedures.

Warnings and Precautions:

- For *In vitro* diagnostic use.
- This control contains human source material and as with all human based specimens, should be treated as potentially infectious and handled with proper laboratory safety procedures to minimize the risk of transmission of infectious disease.
- All components of this kit can be discarded as biohazardous waste according to the local guidelines.
- Refer to the product safety data sheet for risk and safety phrases and disposal information. The product safety data sheet is available via our website, lumiraDx.com.
- Requirements of the appropriate handling or acclimation body should be incorporated into your quality control program.
- Exercise the normal precautions required for handling all laboratory reagents.

Storage and stability:

- Store controls between 2°C and 8°C (36 - 46°F).
- Unopened controls that are stored between 2°C and 8°C (36 - 46°F) can be used until the expiration date. Do not use kits or components beyond the expiration date given on the label.
- Once opened, the vial has a 30 day expiry.
- Open the Control Vials only when you are performing tests.
- Recap and store the Control Vials in their original container at 2°C and 8°C (36 - 46°F) after use.

Carton Contents:

- 2 x 0.5mL vial SARS-CoV-2 & RSV Positive Control
- 2 x 0.5mL vial SARS-CoV-2 & RSV Negative Control
- 40 Transfer pipettes (20µL)
- LumiraDx SARS-CoV-2 & RSV Quality Control Pack Insert

Materials required but not provided with the Control Carton:

- LumiraDx Instrument
- LumiraDx SARS-CoV-2 & RSV Test Strips
- LumiraDx Connect - if connectivity required (refer to LumiraDx Connect User Manual)

Getting ready to test: You will need the LumiraDx Instrument and the following supplies:

- LumiraDx SARS-CoV-2 & RSV Test Strip(s)
- LumiraDx SARS-CoV-2 & RSV Positive or Negative Quality Controls
- Transfer pipette

Handling the LumiraDx SARS-CoV-2 & RSV Test Strips: To ensure that you are using the SARS-CoV-2 & RSV Test and the Instrument correctly, read the appropriate SARS-CoV-2 & RSV Test Strip Product Insert and Platform User Manual.

Preparing the Quality Controls: The liquid controls are supplied ready to use.

Handing the LumiraDx SARS-CoV-2 & RSV Test Strips: To ensure that you are using the SARS-CoV-2 & RSV Test and the Instrument correctly, read the appropriate SARS-CoV-2 & RSV Test Strip Product Insert and Platform User Manual.

Preparation of the Quality Controls: The liquid controls are supplied ready to use.

Handing the LumiraDx SARS-CoV-2 & RSV Test Strips: To ensure that you are using the SARS-CoV-2 & RSV Test and the Instrument correctly, read the appropriate SARS-CoV-2 & RSV Test Strip Product Insert and Platform User Manual.

Preparation of the Quality Controls: The liquid controls are supplied ready to use.

Handing the LumiraDx SARS-CoV-2 & RSV Test Strips: To ensure that you are using the SARS-CoV-2 & RSV Test and the Instrument correctly, read the appropriate SARS-CoV-2 & RSV Test Strip Product Insert and Platform User Manual.

LumiraDx customer services: For product inquiries please contact LumiraDx Customer Services by email: customerservices@lumiraDx.com Any adverse results experienced with the use of this product, and/or quality problems should also be reported to LumiraDx Customer Services by email: customerservices@lumiraDx.com or at lumiraDx.com

For return policy: If there is a problem with the LumiraDx SARS-CoV-2 & RSV Quality Controls Pack you may be asked to return them. Before returning kits, please obtain a return authorization number from LumiraDx Customer Services. The return authorization number must be on the shipping carton for return. For ordinary returns following purchase, please contact LumiraDx Customer Services for terms and conditions.

Procedure/Performing a Test: Consult the LumiraDx Platform User Manual for instructions on how to analyse a Quality Control sample. Open the foil pouch of the SARS-CoV-2 & RSV Test Strip just before use and insert the Test Strip into the LumiraDx Instrument. Hold the Test Strip by gripping the blue label end with the label facing upward. Do not touch Test Strip Sample Application Area. Do not bend or fold the Test Strip. Do not touch Test Strip contacts. The Instrument will indicate when ready for the sample to be applied.

1. Draw up the Quality Control solution into the Transfer pipette. Squeeze the bulb, drawing up QC material to where the tip of the pipette tapers as indicated by the arrow on the diagram. Take care to fill only the tip of the pipette to collect approximately 20µL of QC material.

2. Apply the Quality Control solution to the already inserted SARS-CoV-2 & RSV Test Strip. Hold the pipette over the Sample Application Area of the Test Strip and dispense the Quality Control solution. The LumiraDx Instrument will indicate sample is detected with an audible alert (if the Instrument sounds are enabled). The screen of the LumiraDx Instrument will request the user to close the door. Dispose of the pipette.

3. Do not open the door while the test is in progress. The touch-screen will indicate test progress.

4. The Result will appear on the touchscreen within 12 minutes of applying the Quality Control solution and starting the test. The results will be displayed as a **PASS** or **FAIL** QC result for SARS-CoV-2 & RSV on the instrument screen.

NOTE: If you need to repeat a test, use a new Test Strip.

Expected results: The Instrument displays the result as Pass or Fail. The result is automatically saved in the memory of the Instrument. The system is working properly, and all handling has been completed correctly when the test results obtained are displayed as PASS.

Symbol glossary:

Symbol **Meaning**

IVD In Vitro Diagnostic Medical Device

REF Catalogue Number

LOT Batch code/lot Number

UK CA UK conformity assessed under the Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002)

CE "CE Mark": This product fulfils the requirements of the European Directive 98/79/EC on *in vitro* diagnostic medical devices.

UDI Indicates a carrier that contains unique device identifier information

CONTROL Indicates a control material that is intended to verify the performance characteristics of the LumiraDx Instrument. An indication of the level of control is shown in the box, e.g. 1 or 2.

CONTROL- Indicates a negative control

CONTROL+ Indicates a positive control

BIO Indicates that there are potential biological risks associated with the medical device

EC REP Indicates the authorized representative in the European Community/ European Union

Limited warranty: LumiraDx SARS-CoV-2 & RSV Quality Controls Pack - As per shell file.

For the applicable warranty period, LumiraDx warrants that each product shall be (i) of good quality and free of material defects, (ii) function in accordance with the material specifications referenced in the pack insert and (iii) approved by the proper governmental agencies required for the sale of products for their intended use (the "limited warranty"). If the product fails to meet the requirements of the limited warranty, then on the customer's sole remedy, LumiraDx shall either repair or replace, at LumiraDx's discretion, the LumiraDx SARS-CoV-2 & RSV Quality Controls Pack. Except for the limited warranty stated in this section, LumiraDx disclaims any and all warranties, express or implied, including but not limited to, any warranty of merchantability, fitness for a particular purpose and non-fringement regarding the product. LumiraDx's maximum liability with any customer claim shall not exceed the net product price paid by the customer. Neither party shall be liable to the other party for special, incidental or consequential damages, including without limitation, loss of business, profits, data or revenue, even if a party receives notice in advance that these kinds of damages might result.

The Limited Warranty above shall not apply if the customer has subjected the LumiraDx SARS-CoV-2 & RSV Test Strip and Controls Pack to physical abuse, misuse, abnormal use, use inconsistent with the LumiraDx Platform User Manual, SARS-CoV-2 & RSV Test Product Insert or SARS-CoV-2 & RSV Quality Controls Pack insert, fraud, tampering, unusual physical stress, negligence or accidents. Any warranty claim by Customer pursuant to the Limited Warranty shall be made in writing within the applicable Limited Warranty period.

Intellectual property: The LumiraDx Instrument, Test Strips and all provided LumiraDx documentation ("Products") are protected by law. The Intellectual Property of the LumiraDx Products remains at LumiraDx. Details of relevant Intellectual Property regarding our products can be found at lumiraDx.com/IP.

Legal notices: Copyright © 2022 LumiraDx UK and affiliates. All rights reserved. LumiraDx and Flame logo are protected trademarks of LumiraDx International Ltd. Full details of these and other registrations of LumiraDx can be found at lumiraDx.com/IP. All other trademarks are the property of their respective owners.

SPEC-34881 R2 **ART-01682 R2** **Date of Rev 2022-05**

CE CE mark applies to LumiraDx Instrument, Test Strips, Quality Controls, swabs and Connect Hub only.

EC REP LumiraDx AB, Västra Vågen SA, 16961 Soana, Sweden

SPEC-34881 R2 **ART-01682 R2** **Date of Rev 2022-05**

UK CA UK conformity assessed under the Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002)

CE "CE Mark": This product fulfils the requirements of the European Directive 98/79/EC on *in vitro* diagnostic medical devices.

UDI Indicates a carrier that contains unique device identifier information

CONTROL Indicates a control material that is intended to verify the performance characteristics of the LumiraDx Instrument. An indication of the level of control is shown in the box, e.g. 1 or 2.

CONTROL- Indicates a negative control

CONTROL+ Indicates a positive control

BIO Indicates that there are potential biological risks associated with the medical device

EC REP Indicates the authorized representative in the European Community/ European Union

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Mentions légales : Copyright © 2022 LumiraDx UK et sociétés affiliées. Tous droits réservés. LumiraDx et son logo en forme de flamme sont des marques de commerce protégées de LumiraDx International Ltd. Tous les détails sur ces marques de LumiraDx et de ses engagements se trouvent sur le site lumiraDx.com/IP. Tous les autres marques de commerce appartiennent à leurs propriétaires respectifs.
Copyright © 2022 LumiraDx UK e affiliati. Tutti i diritti riservati. LumiraDx è il logo con la fiamma sono marchi protetti di LumiraDx International LTD. Tutti i dettagli di queste e altre registrazioni di LumiraDx sono disponibili sul sito lumiraDx.com/#IP tutti gli altri marchi commerciali sono di proprietà dei rispettivi titolari.
Copyright © 2022 LumiraDx UK y filiales. Todos los derechos reservados. LumiraDx y el logotipo de la fiamma son marcas comerciales protegidas de LumiraDx International LTD. Para obtener información completa sobre estos y otros registros de LumiraDx, visite LumiraDx.com/#IP. Todas las demás marcas comerciales pertenecen a sus respectivos propietarios.
Copyright © 2022 LumiraDx UK y filiales. Todos los derechos reservados. LumiraDx y el logotipo de la fiamma son marcas comerciales protegidas de LumiraDx International LTD. Para obtener información completa sobre estos y otros registros de LumiraDx, visite LumiraDx.com/#IP. Todas las
LumiraDx UK Ltd, Dumyat Business Park, Alton, FK10 2PB, UK. Número de registro: 092906123
EC REP LumiraDx AB, Västra Vägen 5A, 16961 Solna, Svezia