 **LumiraDx™ HbA1c**

For Professional Use Only

SPEC-35864 R6 ART-02670 R6 Date of Rev 2024-02

**IVD**

Product Name	REF	$\Sigma$
LumiraDx HbA1c	L006000101012	12
LumiraDx HbA1c	L006000101024	24
LumiraDx HbA1c	L006000101048	48

<b>ENGLISH</b>	01 - 17
<b>NEDERLANDS</b>	21 - 38
<b>FRANÇAIS</b>	43 - 60
<b>DEUTSCH</b>	65 - 82
<b>ITALIANO</b>	87 - 103
<b>ESPAÑOL</b>	109 - 126

## For Professional Use Only

SPEC-35864 R6 ART-02670 R6 Date of Rev 2024-02

IVD

### LumiraDx HbA1c:

The LumiraDx HbA1c Test Strips (hereafter referred to as Test Strips) are to be used with the LumiraDx Platform. The LumiraDx Platform is a point of care system for professional use which is used for *in vitro* diagnostic tests. It comprises a portable LumiraDx Instrument, lysis device and a LumiraDx Test Strip for the required test. This test is for **HEALTHCARE PROFESSIONAL USE ONLY** and allows users to perform tests and to view results quickly on the Instrument touch-screen.

### Intended use:

The LumiraDx HbA1c test is an *in vitro* diagnostic test for the quantitative determination of hemoglobin A1c (IFCC mmol/mol and NGSP %) in human capillary and venous whole blood samples (EDTA). The LumiraDx HbA1c Test Strips are intended for use with the LumiraDx Instrument. It is an automated *in vitro* diagnostic test for near-patient testing. HbA1c is used for the monitoring of long-term glycaemic control in individuals with diabetes mellitus, and as an aid in screening and identifying patients who may be at risk for developing diabetes. The LumiraDx HbA1c test is for Professional Use Only. For patients  $\geq 2$  years of age.

**Caution:** For *in vitro* diagnostic use.

Before you start testing, if you are new to the LumiraDx Instrument and LumiraDx Platform, you must read the LumiraDx Platform User Manual, the LumiraDx HbA1c Quick Reference Instructions available online and this entire product insert. In addition, please watch the LumiraDx Platform Training Video available at [lumiradx.com](http://lumiradx.com).

### Summary and explanation of the test:

HbA1c is the stable glucose adduct to the N-terminal group of the beta-chain of HbA0. The measurement of HbA1c in human blood is most important for the long-term control of the glycaemic state in diabetic patients<sup>1</sup>.

Diabetes is characterized by chronic hyperglycemia and causes long-term complications such as retinopathy, neuropathy, and nephropathy. It generally accelerates macrovascular and microvascular changes and has been characterised as a global epidemic<sup>2</sup>.

Since the landmark Diabetes Control & Complications Trial (DCCT) & UK Prospective Diabetes Study (UKPDS) clinical trials clearly demonstrated the relationship between glycaemic control, HbA1c, and diabetic complications, HbA1c has been widely used for the routine monitoring of long-term glycaemic status in both type I and type II diabetic patients<sup>2,3,4</sup>.

A HbA1c measurement reflects time-averaged blood glucose during the previous 2-3 months and is routinely used as the gold standard for long-term follow-up and monitoring of glycaemic control<sup>5</sup>.

## Principle of the assay:

The LumiraDx HbA1c test is a rapid microfluidic immunofluorescence assay for use with the LumiraDx Instrument for the quantitative measurement of HbA1c in human capillary and venous whole blood. The test procedure involves the lysis of capillary or venous whole blood, followed by the addition of the hemolysate to the Sample Application Area of the Test Strip. The Instrument is programmed to perform the analysis when the sample has reacted with the reagents. The analysis is based on the amount of fluorescence the Instrument detects within the measurement area of the Test Strip. The concentration of the analyte in the sample is proportional to the fluorescence detected. The results are displayed on the Instrument touch-screen approximately 5 minutes from sample addition.

## Materials provided:

- LumiraDx HbA1c Test Strips packed individually in sealed desiccant foil pouches
- Lysis Devices – lysis buffer in dropper tube with cap, plus transfer device (single use)
- LumiraDx HbA1c Product Insert
- RFID (Radio frequency ID) Tag held inside the Test Strip carton

## Materials required but not provided with the Test Strip carton:

- LumiraDx Instrument
- LumiraDx HbA1c Quality Controls (as required to meet local and organisational compliance)
- Standard blood collection equipment (i.e. high flow lancets, venepuncture, biohazard disposal)
- Calibrated laboratory pipette (able to accurately dispense 15µL of blood)
- LumiraDx Connect if connectivity required (refer to LumiraDx Connect User Manual)

## Reagents:

The LumiraDx HbA1c test contains mouse monoclonal anti-Hb and anti-HbA1c antibodies, magnetic nanoparticles and a fluorescent latex label. The on-board control (OBC) reagent contains fluorescent latex pre-bound to magnetic nanoparticles and is included to confirm test functionality.

## Warnings and precautions:

- For *in vitro* diagnostic use only.
- Do not open the Test Strip until ready for immediate use.
- Discard and do not use any damaged or dropped Test Strips or other materials.
- Inadequate or inappropriate sample collection, storage, and transport can result in incorrect results.
- The test cannot be visually interpreted; the LumiraDx Instrument must be used to generate results.
- Do not use the kit components beyond the expiration date.
- Do not reuse any kit components.
- Samples must be processed as indicated in the Performing a Test section of this Product Insert. Failure to follow the instructions for use can result in inaccurate results.
- All components of this kit should be discarded as Biohazard waste according to local regulations and procedures.
- Refer to the product safety data sheet for risk and safety phrases and disposal information. The product safety data sheet is available via our website [lumiradx.com](http://lumiradx.com).

- Exercise the normal precautions required for handling all laboratory reagents. Wear protective clothing such as laboratory coats, disposable gloves, and eye protection when samples are collected and evaluated.

Proper laboratory safety techniques should be followed at all times when working with patient samples. Patient samples, used Test Strips, used Lysis Devices and used blood collection equipment may be potentially infectious. Proper handling and disposal methods should be established by the laboratory in accordance with local regulations and procedures.

### **Storing the Test Strips (and Lysis Devices):**

Store the Test Strips and Lysis Devices in their original carton. You can store the Test Strips and Lysis Devices at a temperature between 2°C and 30°C (36°F and 86°F). Avoid freezing or storing in any area that could exceed 30°C. When stored properly, the Test Strips and Lysis Devices can be used until the expiration date printed on the Test Strip foil pouch and the Test Strip carton. Discard the Test Strips and Lysis Devices if they are past the expiration date.

### **Handling the Test Strips:**

When you are ready to perform a test, open the Test Strip carton, take out a Test Strip, and remove it from the foil pouch. Hold the Test Strip by gripping the blue label end with the label facing upwards. Do not touch the Test Strip Sample Application Area. Do not bend or fold the Test Strip. Do not touch the Test Strip contacts. After removing the Test Strip from the foil pouch, it should be used immediately. Do not use the Test Strip if there are any visible signs of damage to the foil pouch such as tears or holes.

### **Sample material:**

The following samples can be used with the LumiraDx HbA1c Test Strip:

- Capillary whole blood from a fingerstick sample, delivered via Lysis Device
- Anticoagulated venous whole blood (EDTA), delivered via Lysis Device
- LumiraDx HbA1c Quality Controls, delivered via transfer pipette

### **The test device contains:**

- Mouse monoclonal antibodies
- Fluorescent particles
- Magnetic particles
- Buffer and Stabilising Agents

### **Sample collection and preparation for analysis:**

When collecting any type of sample, follow universal blood collection precautions and guidelines according to your organization. For sample collection of venous whole blood, follow the sample tube manufacturer's recommended procedure.

- Capillary blood samples are stable for up to 5 hours once placed in the lysis buffer solution
- Venous blood should be tested within 24 hours of sample collection or:
  - Refrigerated and tested within 6 days
  - Once placed in the lysis buffer solution, samples are stable for up to 5 hours

### **Preparing the Instrument to perform a test:**

Power on the Instrument by pressing the power button at the rear of the Instrument. You will hear the Instrument powering on, and the display will be a blank black screen for several seconds before starting up. If the screen is just dimmed tap the touch-screen to wake up the Instrument.

Refer to the section on Performing a Test in this Product Insert for information on how to test a Patient sample. The LumiraDx Quick Reference Instructions (QRI) provide an illustrated step-by-step procedure on how to run a Test. Operate the LumiraDx Platform with the HbA1c test at room temperature between 15°C and 30°C (59°F and 86°F) and 10% - 90% relative humidity.

The Instrument will prompt to install the Lot Calibration File when inserting a new Test Strip Lot. Once installed, the Instrument will have all the information required to process the test, and any future tests from the same Lot of Test Strips.

### Lot Calibration File installation:

Lot Calibration Files are required to provide the Instrument with the information needed to perform diagnostic tests. This only needs to be completed once for each Test Strip Lot. The Instrument will prompt to install the Lot Calibration File when inserting a new Test Strip Lot.

#### RFID strip code reader


Locate ((•)) symbol on Instrument.

#### Installation

Place and hold the back of the Test Strip Carton ((•)) symbol to install.

The Instrument will sound and a confirmation message will be displayed.


When indicated by the touchscreen, open the foil pouch just before use and insert the LumiraDx Test Strip into the LumiraDx Instrument. The Instrument will indicate when it is ready for the sample to be applied.

The LumiraDx HbA1c test results should be evaluated by a Healthcare Professional in the context of all available clinical and laboratory data.



### Testing from a fresh capillary fingerstick sample using the Lysis Device:

The steps that follow apply to collecting a capillary blood sample from a finger stick and applying the sample via the Lysis Device. Only auto-disabling, single use, high flow lancing devices may be used to collect capillary blood.

1. Increasing the blood flow in the finger will help to get a good drop of blood. Before lancing the finger, the following techniques can be used until the fingertip has increased colour:
  - Ask the patient to rinse their hands with warm water.
  - Ask the patient to hold his or her arm straight down at their side
  - Massage the finger from its base, and if required, immediately after lancing, very gently squeeze the finger from its base to encourage blood flow.
2. **Use a high flow lancet** on the selected finger to obtain a large, beaded blood sample, roughly 5mm wide (15µL). Note: It is recommended to use a lower numbered gauge (G) high flow lancet to extract the appropriate blood drop size.

3. Once a large, beaded drop has formed, lightly press and hold the loop of the Transfer Device against the blood drop until the loop is full (see images 1 - 3 below for guidance). Once the loop is full, the blood drop will be held in place by surface tension. Note: please ensure to completely fill the Transfer Device before transferring to the Buffer Tube.


**Note:** If the blood drop runs down the finger, the sample won't collect properly in the loop of the transfer device. Refer to the below images for examples of an ideal blood drop.


**Image 1.**  
Incorrect - the blood drop is running down the finger




**Image 2.**  
Correct - the blood drop is beaded and about 5mm in diameter

4. Place the transfer loop containing the sample into the lysis buffer tube and roll the loop at least 10 times between the thumb and index finger to lyse the blood.




5. Direct the opening of the buffer tube away from the face and snap the transfer loop stem at the break point by pressing it against the inside of the buffer tube and dispose of the upper part of the stem. The loop and remaining stem will stay in the tube with the hemolysate sample.
6. Place the dropper lid on the tube.  
**Note:** The sample is stable at room temperature in the buffer solution for up to 5 hours. Place in a tube holder until ready to run the test.
7. **Apply the hemolysate sample** within 5 hours from sample lysis by positioning the dropper of the Lysis Device and squeezing to generate a single, large drop (20 $\mu$ L) over the Sample Application Area of the inserted Test Strip. Allow the hemolysate drop to touch the Sample Application Area of the Test Strip.
8. Sample will then be drawn by capillary action into the Test Strip. When the sample is detected, the Instrument will sound (if sounds are enabled) and a confirmation message will be displayed. The touch-screen of the LumiraDx Instrument will request the user to close the door.
9. **Do not add more sample.** Do not open the door while the test is in progress. The touch-screen will indicate test progress.
10. **The result** will appear on the Instrument touch-screen within approximately 5 minutes of applying the sample and starting the test.
11. **Dispose** of the Lysis Device, lancet and Test Strip in the appropriate clinical waste.
12. **Clean** the patient's finger with a clean tissue and apply slight pressure.
13. If you need to retest, use a new Test Strip, lancet and Lysis Device, and a different finger.



### Testing from venous whole blood sample using the Lysis Device:

The steps that follow apply to collecting a venous whole blood sample and applying the sample via the Lysis Device. A calibrated laboratory pipette should be used to transfer the blood sample from the venous whole blood tube to the lysis buffer tube.

1. Only use EDTA anticoagulated venous whole blood.
2. Mix the sample well before testing.
3. Carefully remove the lid of the venous whole blood tube.
4. Using a calibrated laboratory pipette (capable of accurately delivering 15 $\mu$ L of blood) or equivalent device, aspirate 15 $\mu$ L of blood into the tip of the pipette. Be sure to check for excess blood on the outside of the pipette tip. Wipe away any excess blood if necessary.

5. Dispense 15 $\mu$ L of blood from the venous whole blood tube into the lysis buffer tube.
6. Dispose of the pipette tip in the appropriate clinical waste.
7. Place the transfer loop into the lysis buffer tube and roll the loop at least 10 times between the thumb and index finger to lyse the blood.
8. Direct the opening of the buffer tube away from the face and snap the transfer loop stem at the break point by pressing it against the inside of the buffer tube and dispose of the upper part of the stem. The loop and remaining stem will stay in the tube with the hemolysate sample.


9. Place the dropper lid on the tube.  
**Note:** The sample is stable at room temperature in the buffer solution for up to 5 hours. Place in a tube holder until ready to run the test.
10. **Apply the hemolysate sample** within 5 hours from sample lysis by positioning the dropper of the Lysis Device and squeezing to generate a single, large drop (20 $\mu$ L) over the Sample Application Area of the inserted Test Strip. Allow the hemolysate drop to touch the Sample Application Area of the Test Strip.

11. Sample will then be drawn by capillary action into the Test Strip. When the sample is detected, the Instrument will sound (if sounds are enabled) and a confirmation message will be displayed. The touch-screen of the LumiraDx Instrument will request the user to close the door.
12. **Do not add more sample.** Do not open the door while the test is in progress. The touch-screen will indicate test progress.
13. **The result** will appear on the Instrument touch-screen within approximately 5 minutes of applying the sample and starting the test.
14. **Dispose** of the Lysis Device, and Test Strip in the appropriate clinical waste and clean all surfaces as per local guidelines.
15. If you need to retest, use a new Test Strip, and Lysis Device.

### Results interpretation:

The results will be displayed on the Instrument screen –  
**example of results screen display:**

## Expected Values:

Interpretation of results may vary slightly by region or by guideline used. The below was derived from the American Diabetes Association (ADA) guidelines<sup>6</sup>:

- Non-diabetic: <39 mmol/mol (5.7% HbA1c)
- Pre-diabetic: 39 - 47 mmol/mol (5.7 - 6.4% HbA1c)
- Diabetic:  $\geq$  48 mmol/mol ( $\geq$ 6.5% HbA1c)

Each laboratory should establish a reference range that is representative of the patient population to be evaluated. Additionally, each laboratory should consider the current practice in the evaluation of patients experiencing symptoms at each institution.


The LumiraDx HbA1c test results should be evaluated by a healthcare professional in the context of all available clinical and laboratory data.

## Invalid test results:

If an issue occurs, a message will be displayed on the Instrument touch-screen. Alert messages include useful information and are highlighted by an orange banner. Error messages also include a  symbol. All messages will contain a description of the Instrument status or error and an instruction. Error messages contain an identifying code that may be used for further troubleshooting purposes. Refer to the LumiraDx Platform User Manual if an error message is displayed on the LumiraDx Instrument touch-screen and contact LumiraDx Customer Services.

### Example of an error screen:

If the On Board Control (OBC) fails, an error message will be shown and no test result will be returned. Follow the on screen instructions to dispose of the Test Strip and start a new test. If the problem persists, contact Customer Services.

## Testing patient samples procedural notes:

- Refrigerated whole blood samples must be allowed to reach room temperature and be mixed thoroughly before testing.
- Before use, mix whole blood venous samples by gently inverting the tube several times.

## Built-in controls:

The Instrument reads the 2D barcode on each Test Strip and can identify if the Test Strip has exceeded the expiry date for use, and if the Test Strip Lot Calibration file has not yet been loaded at which point it will request it.

The LumiraDx Instrument and LumiraDx HbA1c Test Strips have several quality control functions integrated to ensure validity of each test run. These checks ensure that the volume of sample added is sufficient and the assay sequence of the Test Strip is as expected. The checks also ensure that the Test Strip has not been damaged or used previously. If these checks are not verified, the test run will be rejected and an error message displayed on the Instrument touchscreen.

The LumiraDx Instrument ensures the quality of test results obtained through the following features:

- Automated checks of the correct functioning of the Instrument at power on and during operation. This includes electrical component operation, heater operation, battery charge state, mechanical actuators and sensors and optical system performance.
- Monitoring of Test Strip performance and controls during test runtime.
- Ability to perform Quality Control Tests using LumiraDx Quality Control solutions to meet regulatory compliance requirements.

## Standardisation:

The calibration of the LumiraDx HbA1c test is traceable to the International Federation of Clinical Chemistry and Laboratory Medicine (IFCC) primary reference method for the measurement of HbA1c.

## External Quality Controls:

External liquid Quality Controls for LumiraDx HbA1c are available from LumiraDx and may be used to demonstrate that the Test is functioning properly by demonstrating the expected Quality Control results and correct test performance by the operator. External Quality Control requirements should be established in accordance with local and organizational compliance. It is recommended that external control testing be performed with each new operator and before using a new lot or shipment of the LumiraDx HbA1c. Refer to the LumiraDx HbA1c Quality Controls pack insert available at [lumiradx.com](http://lumiradx.com) for detailed instructions. LumiraDx HbA1c Quality Controls are purchased separately.

If the LumiraDx HbA1c Quality Controls do not perform as expected, repeat the QC Test and if the problems persist, do not report patient results and contact LumiraDx Customer Services.

## Cleaning and disinfection:

It is recommended to disinfect the Instrument after each patient sample, or if contamination is suspected. Excessive liquid may damage the Instrument. It is important for the protection of the Instrument that exposure to excess moisture is prevented. All disinfection cloths and/or wipes should only be slightly damp, with any excess liquid being manually removed from the cloth before use. Alcohol wipes alone are not sufficient to disinfect the Instrument for blood-based samples, due to the potential presence of bloodborne pathogens

1. Using a LumiraDx recommended disinfecting material, wipe the external surfaces of the Instrument while taking care to avoid the door hinges, Test Strip inlet, power cord, and USB port.
2. Allow the disinfectant at least 5 minutes contact time with the Instrument before testing the next sample.
3. Dispose of disinfectant materials in accordance with local biohazardous waste disposal procedures.

To clean the Instrument wipe the external surfaces with a soft, slightly damp cloth when it appears visibly dirty.

For more information, or for the full procedure on cleaning and disinfection, please refer to the Technical Bulletin Platform Disinfection Procedure at [lumiradx.com](http://lumiradx.com).

## Limitations:

- The LumiraDx HbA1c test uses fresh capillary whole blood (15 $\mu$ L) and venous blood (15 $\mu$ L) which is collected and applied via a Lysis Device. The resultant hemolysate sample size must be a minimum of 20 $\mu$ L in volume. Low sample volume will cause an error message. Never add more sample to the Test Strip after the test has begun.
- Use the Test Strip only once and then dispose of it appropriately in clinical waste.
- There is the possibility that factors such as technical or procedural errors, as well as additional substances in blood samples that are not listed below, may interfere with the test and cause erroneous results.
- Blood sample types, draw methods or anticoagulants different from those described in this pack insert have not been evaluated.
- As with any assay employing mouse antibodies, the possibility exists for interference by human anti-mouse antibodies (HAMA) in the sample. The test has been formulated to minimize this interference; however, samples from patients who have been routinely exposed to animal serum products may contain heterophile antibodies which may cause erroneous results.

- The test has been formulated to minimise interference from Rheumatoid Factors (RF), however, due to heterogeneity of RF, specimens from patients with highly elevated RF may cause erroneous results.
- Any unusual result must always be followed up to identify the potential cause.
- Marked discordance between HbA1c and plasma glucose levels should be followed up to identify the potential cause<sup>6</sup>.
- In conditions associated with an altered relationship between HbA1c and glycemia and/or decreased red blood cell survival, the LumiraDx HbA1c test should not be used. Some examples include disease-state/homozygous hemoglobinopathies (including sickle cell disease), hemolytic anemia, iron deficiency anemia, polycythemia, pregnancy (second and third trimesters and the postpartum period), glucose-6-phosphate dehydrogenase deficiency, malignancies, severe chronic hepatic and renal diseases, HIV, hemodialysis, recent blood loss, recent transfusion, or erythropoietin therapy.
- Specimens containing HbF (>7.2%), including patients with Hereditary Persistence of Fetal Hemoglobin (HPFH), may result in lower than expected HbA1c values and should not be reported.
- Results that do not match the clinical symptoms should be repeated to rule out a procedural error.
- The assay has not been validated for individuals younger than 2 years of age.
- When performing a new test or repeating a patient test, always use a new lancet to obtain a fresh drop of blood from a different finger and use a new Test Strip and Lysis Device.
- The LumiraDx HbA1c test is only validated to analyze fresh capillary and refrigerated venous whole blood samples. Please refer to the sample collection and preparation for analysis section for further information. Frozen samples are not a validated sample matrix for this test.

Unusual Results: If the LumiraDx Instrument displays an error message, refer to the Troubleshooting section of the LumiraDx Platform User Manual. If the LumiraDx Instrument displays an unexpected test result (other than an error message), check this Limitations section.

## Results:

The LumiraDx HbA1c test measures HbA1c concentration via measurement of an optical signal generated when the fluorescent immunoassay (FIA) reagents deposited on the Test Strip are resuspended in the patient sample. The measured optical signal is proportional to the HbA1c concentration. The optical signal is then converted to HbA1c concentration via use of a calibration curve, which is established per lot of Test Strips during the calibration process.

Each HbA1c result is reported on screen in any 1 of the following user-selectable options:

1. Both mmol/mol (IFCC), and %HbA1c (NGSP) units of measurement
2. mmol/mol (IFCC) units of measurement only
3. %HbA1c (NGSP) units of measurement only

The system default unit of measurement is both mmol/mol (IFCC) and %HbA1c (NGSP) units of measurement. Units of measurement are configurable via the Settings Menu. Please refer to the Platform User Manual for more information.

The IFCC-NGSP master equation is used for conversion between units of measurement:  $\text{HbA1c (\%)} = \text{IFCC} \times 0.09148 + 2.152^4$ .

## Performance characteristics:

### Measuring range:

The LumiraDx HbA1c test used with the LumiraDx Instrument has a reportable range of 20 - 130 mmol/mol HbA1c (4.0 - 14.0% HbA1c). HbA1c <20 mmol/mol HbA1c (<4.0% HbA1c) is displayed if the HbA1c concentration is less than 20 mmol/mol HbA1c (4.0% HbA1c). HbA1c >130 mmol/mol HbA1c (>14.0% HbA1c) is displayed if the HbA1c concentration is greater than 130 mmol/mol HbA1c (14.0% HbA1c).

The hemoglobin (Hb) measuring range for the LumiraDx HbA1c Test is 6-21 g/dL.

**Note:** hemoglobin (Hb) is not reported as part of the LumiraDx HbA1c Test.

#### Linearity:

Linearity was established according to a protocol based on CLSI EP06 ED2<sup>7</sup> in venous whole blood (EDTA). High HbA1c concentration samples were sourced from diabetic clinics. Linearity series were then prepared by mixing the high concentration samples with blood depleted of analyte. The results obtained, confirm linearity across the measuring range of 20 - 130 mmol/mol (4.0 - 14.0%).

#### Hook effect:

No hook effect is observed with the LumiraDx HbA1c test at HbA1c concentrations up to 151 mmol/mol HbA1c (16.0% HbA1c).

#### Precision:

A precision study was carried out in venous whole blood (EDTA) on a protocol based on CLSI EP05-A3<sup>8</sup>. The study was carried out at 2 concentrations of HbA1c, each was tested in 1 run of 5 replicates per day, for 5 days across 3 sites. The results of the precision study are summarised below:

HbA1c Concentration mmol/mol HbA1c (% HbA1c)	Within Day Precision (%CV)	Between Day Precision (%CV)	Between Site Precision (%CV)	Total Precision (%CV)	n
Level 1 (6.5-7.5%)	2.9	2.0	1.9	4.0	70
Level 2 (8.5-9.5%)	2.3	1.7	1.6	3.4	75

Capillary & venous whole blood precision was also determined using duplicate samples from patients tested on 1 Test Strip lot at multiple sites across a HbA1c range of 31.4 - 128.7 mmol/mol HbA1c (5.1 - 12.8% HbA1c). The following results represent the mean paired rep %CV (mprCV) for each sample type tested:

Sample Type	NGSP (%)			IFCC (mmol/mol)		
	N	Range	mprCV	N	Range	mprCV
Capillary Blood	62	5.1 - 12.1	2.50	62	32.1 - 117.5	3.40
Venous Blood	64	5.1 - 12.8	2.22	63	31.4 - 128.7	3.12

#### Method comparison:

The method comparison was performed using 2 Test Strip lots with fresh whole blood samples sourced from patients presenting to outpatient clinics or research facilities with suspected or confirmed diabetes, or from normal healthy donors.

Each sample tested on the LumiraDx Platform was compared to the same fresh whole blood sample tested on the Tosoh G8 (Tosoh Automated Glycohemoglobin Analyzer HLC-723G8). The data was analysed by Passing Bablok regression. The analyses are summarised below:

LumiraDx Test Strip Lot	n	HbA1c range mmol/mol (% HbA1c)	Slope	Intercept	r
NGSP all lots	394	(4.8 - 13.5)	1.06	-0.26	0.98
IFCC all lots	389	28.0 - 129.0	1.03	-0.97	0.97
NGSP (Lot 1)	196	(4.8 - 13.4)	1.05	-0.30	0.98
NGSP (Lot 2)	198	(5.1 - 13.5)	1.06	-0.24	0.98
IFCC (Lot 1)	196	28.0 - 126.0	1.05	-1.69	0.97
IFCC (Lot 2)	193	33.0 - 129.0	1.00	1.00	0.97

#### Matrix equivalency:

A study was carried out with 57 subjects presenting to outpatient clinics or research facilities with suspected or confirmed diabetes. Samples of capillary fingerstick blood (via Lysis Device application) and paired whole blood (EDTA via Lysis Device application) were collected and tested. The data were analysed by Passing Bablok regression. The analyses are summarised below:

Reference	N	Range	Slope	Intercept	r
NGSP	57	5.1 - 12.4%	0.98	0.19	0.98
IFCC	57	32.1 - 128.7 mmol/mol	0.99	0.54	0.98

#### Interference:

Testing was performed according to a protocol based on CLSI EP07 ED3<sup>9</sup>. Testing was carried out using EDTA whole blood samples, where possible, at 2 concentrations of HbA1c: 48 - 58 mmol/mol (6.5 - 7.5%) and 69 - 80 mmol/mol (8.5 - 9.5%) spiked with interfering substances. The following interferences showed no significant effect on HbA1c Test results. Significant interference was defined as  $\geq \pm 7\%$  mathematical bias from a reference value for samples  $\leq 68$  mmol/mol (8.4% HbA1c), and  $\geq \pm 10\%$  mathematical bias from a reference value for samples  $\geq 69$  mmol/mol (8.5% HbA1c).

#### Exogenous (test concentration):

Amlodipine (0.0075 mg/dL), Amoxicillin (5.4 mg/dL), Apixaban (0.0315 mg/dL), Ascorbic acid (5.25 mg/dL), Aspirin (3 mg/dL), Atenolol (0.9 mg/dL), Atorvastatin (0.075 mg/dL), Biotin (0.351 mg/dL), Bisoprolol (0.0258 mg/dL), Caffeine (10.8 mg/dL), Calciferol (0.015 mg/dL), Cetirizine (0.435 mg/dL), Citalopram (0.543 mg/dL), Clarithromycin (0.72 mg/dL), Clopidogrel (18 mg/dL), Doxycycline (1.08 mg/dL), EDTA (0.099 mg/dL), Empagliflozin (1.5 mg/dL), Fluconazole (2.55 mg/dL), Folic acid (0.006 mg/dL), Furosemide (1.59 mg/dL), Gliclazide (1.92 mg/dL), Heparin Sodium (1.76 mg/dL), Ibuprofen (21.9 mg/dL), Levothyroxine (0.0429 mg/dL), Linagliptin (0.3 mg/dL), Metformin (5.76 mg/dL), Metronidazole (12.3 mg/dL), Nicotine (0.097 mg/dL), Omeprazole (0.84 mg/dL), Paracetamol (15.6 mg/dL), Prochlorperazine (0.345 mg/dL), Ramipril (0.0156 mg/dL), Salbutamol (0.0045 mg/dL), Sertraline (0.093 mg/dL), Simvastatin (0.168 mg/dL), Warfarin (7.5 mg/dL), Insulin aspart (0.0039 mg/dL), Insulin glargine (0.00404 mg/dL).

**Endogenous (test concentration):**

Bilirubin (conjugated) (74.06 mg/dL), Bilirubin (unconjugated) (16.16 mg/dL), Glucose (1000 mg/dL), Total Protein (16.58 g/dL), Hemolysis (1480 mg/dL), and Lipemia (Triglycerides) (593.4 mg/dL).

It is possible that other substances and/or factors not listed above may interfere with the test and cause inaccurate results.

**Hemoglobin variants and cross-reactants:**

For HbA0, HbA1a, HbA1b, Acetylated Hemoglobin, Carbamylated Hemoglobin, Labile Hemoglobin & Glycated Albumin, a study was performed based on guidance from CLSI EP07 ED3<sup>o</sup>. Specificity was assessed by comparing samples containing abnormal hemoglobins and cross-reactants to a comparator assay. Significant interference was defined as  $\geq \pm 6\%$  from a reference value.

For HbA2, HbC, HbD, HbE, HbF & HbS, a study was performed using fresh whole blood samples containing common hemoglobin variants. The results were compared to a method that has been confirmed free of interference. Significant interference was defined as  $\geq \pm 6\%$  mean relative deviation from a reference value.

Hb Species	Highest level tested with no significant interference
HbA0	93.6%
HbA1a	1.3%
HbA1b	2.0%
HbA2	6.7%
HbC	44%
HbD	40%
HbE	35%
HbF	7.2%
HbS	40%
Acetylated Hemoglobin	2.76 mg/mL
Carbamylated Hemoglobin	13.8 mg/mL
Labile Hemoglobin	11.4 mg/mL
Glycated Albumin	3.85 mg/mL

It is possible that other substances and/or factors not listed above may interfere with the test and cause inaccurate results.

## Symbols glossary

Symbol	Meaning
	Temperature limitation
	Manufacturer
	Importer
	Distributor
<b>IVD</b>	<i>In vitro</i> diagnostic medical device
<b>REF</b>	Catalogue Number
<b>LOT</b>	Lot Number
	Use-by Date - indicates the date after which the unopened IVD/ Quality Control Material cannot be used
<b>CE</b>	"CE Mark". This product fulfils the requirements of the European Directive 98/79/EC on <i>in vitro</i> diagnostic medical devices.
<b>UK CA</b>	UK conformity assessed under the Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002)
	Refer to instructions for use
	Do not re-use
	Indicates the presence of the Radio Frequency Identification (RFID) reader/tag.
<b>EC REP</b>	Authorized Representative in the European Union
	Total number of IVD tests that can be performed with the IVD medical device.
	For near patient testing
<b>UDI</b>	Indicates a carrier that contains Unique Device Identifier information

## References:

1. Jeppsson JO et al. Approved IFCC Reference Method for the Measurement of HbA1c in Human Blood. Clin Chem Lab Med. 2002.
2. Weykamp C. HbA1c: A Review of Analytical and Clinical Aspects. Annals of Laboratory Medicine. 2013.
3. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med 1993.
4. King P et. al. The UK Prospective Diabetes Study (UKPDS): clinical and therapeutic implications for type 2 diabetes. J Clin Pharmacol. 1999.
5. Hanas R, John WG. 2010 Consensus Statement on the Worldwide Standardization of the Hemoglobin A1c Measurement. Diabetes Care. 2010.
6. American Diabetes Association. Standards of Medical Care in Diabetes. 2023.
7. CLSI EP06 ED2
8. CLSI EP05 A3
9. CLSI EP07 ED3

Visit [www.CLSI.org](http://www.CLSI.org) for information

## LumiraDx customer services:

For product inquiries please contact LumiraDx Customer Services at [customerservices@lumiradx.com](mailto:customerservices@lumiradx.com) or find telephone contact details at [lumiradx.com](http://lumiradx.com).

If during the use of the device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

## For return policy:

If there is a problem with the LumiraDx HbA1c tests you may be asked to return them. Before returning tests please obtain a return authorization number from LumiraDx Customer Services. This return authorization number must be on the shipping carton for return. For ordinary returns following purchase, please contact LumiraDx Customer Services for terms and conditions.

## Limited warranty:

LumiraDx HbA1c Test Strips – as per shelf life.

Unused strips must be stored according to the required storage conditions as printed in this product insert and they can be used only up to the expiry date printed on the Test Strip pouch and Test Strip box. For the applicable warranty period, LumiraDx warrants that each product shall be (i) of good quality and free of material defects, (ii) function in accordance with the material specifications referenced in the product insert, and (iii) approved by the proper governmental agencies required for the sale of products for their intended use (the "limited warranty"). If the product fails to meet the requirements of the limited warranty, then as customer's sole remedy, LumiraDx shall either repair or replace, at LumiraDx's discretion, the Test Strips. Except for the limited warranty stated in this section, LumiraDx disclaims any and all warranties, express or implied, including but not limited to, any warranty of merchantability, fitness for a particular purpose and non-infringement regarding the product. LumiraDx's maximum liability with any customer claim shall not exceed the net product price paid by the customer. Neither party shall be liable to the other party for special, incidental or consequential damages, including, without limitation, loss of business, profits, data or revenue, even if a party receives notice in advance that these kinds of damages might result. The Limited Warranty above shall not apply if the customer has subjected the LumiraDx HbA1c Test Strips to physical abuse, misuse, abnormal use, use inconsistent with the LumiraDx Platform User Manual or Product Insert, fraud, tampering, unusual physical stress, negligence or accidents. Any warranty claim by Customer pursuant to the Limited Warranty shall be made in writing within the applicable Limited Warranty period.

### Intellectual property:

The LumiraDx Instrument, Test Strips and all provided LumiraDx documentation ("Products") are protected by law. The Intellectual Property of the LumiraDx Products remains at LumiraDx. Details of relevant Intellectual Property regarding our products can be found at [lumiradx.com/IP](http://lumiradx.com/IP). Test Strips contained herein include reagents provided under license right from TriLink BioTechnologies LLC.

### Legal notices:

Copyright © 2024 LumiraDx UK and affiliates. All rights reserved. LumiraDx and Flame logo are protected trademarks of LumiraDx International LTD. Full details of these and other registrations of LumiraDx can be found at [lumiradx.com/IP](http://lumiradx.com/IP). All other trademarks are the property of their respective owners.

### Manufacturer information:

**LumiraDx UK Ltd**  
Dumyat Business Park, Alloa  
FK10 2PB, UK  
**Registration number:** 09206123

CE mark applies to LumiraDx Instrument, Test Kit, Quality Controls and Connect Hub only.

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1. Jeppsson JO et al. Approved IFCC Reference Method for the Measurement of HbA1c in Human Blood. Clin Chem Lab Med. 2002.
2. Weykamp C. HbA1c: A Review of Analytical and Clinical Aspects. Annals of Laboratory Medicine. 2013.
3. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med 1993.
4. King P et. al. The UK Prospective Diabetes Study (UKPDS): clinical and therapeutic implications for type 2 diabetes. J Clin Pharmacol. 1999.
5. Hanas R, John WG. 2010 Consensus Statement on the Worldwide Standardization of the Hemoglobin A1c Measurement. Diabetes Care. 2010.
6. American Diabetes Association. Standards of Medical Care in Diabetes. 2023.
7. CLSI EP06 ED2
8. CLSI EP05-A3
9. CLSI EP07 ED3

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Copyright © 2024 LumiraDx UK en gelieerde ondernemingen. Alle rechten voorbehouden. LumiraDx en het vlamlogo zijn beschermde handelsmerken van LumiraDx International LTD. Volledige gegevens over deze en andere gedeponeerde handelsmerken van LumiraDx vindt u op [lumiradx.com/IP](http://lumiradx.com/IP). Alle overige handelsmerken zijn eigendom van de desbetreffende merkhouders.
Copyright © 2024 LumiraDx UK und Tochtergesellschaften. Alle Rechte vorbehalten. LumiraDx und das Flammen-Logo sind geschützte Marken von LumiraDx International LTD. Vollständige Angaben zu diesen und anderen Eintragungen von LumiraDx sind zu finden unter [lumiradx.com/IP](https://lumiradx.com/IP). Alle anderen Marken sind Eigentum des jeweiligen Inhabers.
Copyright © 2024 LumiraDx UK e affiliate. Tutti i diritti riservati. LumiraDx e il logo con la fiamma sono marchi protetti di LumiraDx International LTD. Tutti i dettagli di queste e altre registrazioni di LumiraDx sono disponibili sul sito [lumiradx.com/IP](http://lumiradx.com/IP). Tutti gli altri marchi di fabbrica sono di proprietà dei rispettivi titolari.
Copyright © 2024 LumiraDx UK y filiales. Todos los derechos reservados. LumiraDx y el logotipo de la llama son marcas comerciales protegidas de LumiraDx International LTD. Para obtener información completa sobre estos y otros registros de LumiraDx, visite [lumiradx.com/IP](http://lumiradx.com/IP). Todas las demás marcas comerciales pertenecen a sus respectivos propietarios.
LumiraDx B.V. Looskade 20
LumiraDx AB, Västra Vägen 5A,