



D-Dimer Test Specifications

For *in vitro* diagnostic use.

Intended Use

The LumiraDx D-Dimer test is an *in vitro* diagnostic test for the quantitative determination of D-Dimer in human capillary and venous whole blood and plasma samples (Sodium Citrate). The LumiraDx D-Dimer Test Strips are intended for use with the LumiraDx Instrument. It is an automated *in vitro* diagnostic test for near-patient testing to aid in the assessment and diagnosis of patients with suspected venous thromboembolism (VTE) such as deep vein thrombosis (DVT) and pulmonary embolism (PE).

The test can be used in conjunction with a clinical pre-test probability assessment model to exclude deep vein thrombosis (DVT) and pulmonary embolism (PE) disease in patients suspected of DVT or PE. The LumiraDx D-Dimer test is for Professional Use Only. For patients ≥ 18 years of age.

Test Description

D-Dimer is a degradation product of fibrin, present in the blood after a blood clot is degraded by fibrinolysis. D-Dimer testing is of clinical use when there is a suspicion of VTE and is used alongside clinical scoring systems and additional test methods.

The LumiraDx D-Dimer test is an easy to use, fast microfluidic immunoassay designed to rapidly quantify D-Dimer levels in whole blood and plasma.

The test result is the mean of 3 D-Dimer assays run on the unique multi-channel Test Strip.

The LumiraDx D-Dimer test is the only direct fingerstick D-Dimer assay available at the point of care today*, aiding healthcare professionals to exclude deep vein thrombosis (DVT) and pulmonary embolism (PE) in symptomatic patients with confidence - all in only 6 minutes.

*As stated at time of publication - 24th May 2022

Built-in Quality Controls

LumiraDx Platform Instrument and Test Strips are integrated with several control checks to ensure that the Instrument and Test are functioning correctly for every test run. These checks include:

- Electrical component operation, heater operation, battery charge state, mechanical actuators and sensors and optical system performance
- Test Strip positioning and Test Strip expiry
- Monitoring of the Test Strip performance and controls during test runtime
- The D-Dimer test contains on Onboard Quality Control (OBC) assay
- Sufficient sample volume and hematocrit determination on the Test Strip to ensure patients are within the 20-55% range and results compensated for blood performance

D-Dimer Quality Controls

D-Dimer Quality Controls come in two levels and are available from LumiraDx.

Method Comparison

The method comparison was performed using plasma samples from patients ($n = 327$, range = 60 - 4515 $\mu\text{g/L}$ FEU (0.060-4.515 mg/L FEU). A comparison of 1767 D-Dimer measurements with the LumiraDx D-Dimer test to the VIDAS Exclusion II D-Dimer assay yielded the following statistics: Slope = 1.02, Intercept = 21, $r = 0.92$.

Precision

A precision study was carried out in citrated venous plasma on a protocol based on CLSI EP5-A3. The study was carried out with levels of D-Dimer, each was tested in 2 runs of 2 replicates per day, for twenty days.

| D-Dimer concentration ($\mu\text{g/L}$ FEU) | D-Dimer concentration (mg/L FEU) | Total precision (% CV) | n |
|--|--|------------------------|----|
| 291 | 0.291 | 11.1 | 80 |
| 552 | 0.552 | 9.7 | 80 |
| 1790 | 1.790 | 10.2 | 80 |

Clinical Performance

A prospective clinical study was performed on 585 subjects where fresh samples (capillary blood, venous (blood citrated) and plasma (citrated)) were collected from patients presenting with symptoms of VTE (PE or DVT) ^{1,2}. Subjects also required an assessment with the Wells score and were classed as pre-test probability (PTP) 'Likely' or PTP 'Unlikely'. Those with 'Unlikely' PTP categorization were further analysed using the LumiraDx D-Dimer test with 500 $\mu\text{g/L}$ FEU (0.500 mg/L FEU) D-Dimer as cut-off. The corresponding sensitivity and negative predictive values (NPV) by sample matrix are listed below with corresponding Wilson Score 95% confidence intervals (CI).

| Estimate | Matrix | Patients with Suspected VTE |
|------------------------|------------------|--------------------------------|
| | | Unlikely PTP |
| Sensitivity % (95% CI) | Venous | 100.0% (74.1%-100.0%; n = 378) |
| | Capillary Direct | 100.0% (72.2%-100.0%; n = 377) |
| | Plasma | 100.0% (74.1%-100.0%; n = 406) |
| NPV % (95% CI) | Venous | 100.0% (98.3%-100.0%; n = 378) |
| | Capillary Direct | 100.0% (98.1%-100.0%; n = 377) |
| | Plasma | 100.0% (98.1%-100.0%; n = 406) |

D-Dimer Test Strip Specifications

| | |
|-----------------------|--|
| Clinical cut-off | 500 µg/L FEU (0.500 mg/L FEU) |
| Displayed results | D-Dimer FEU (Fibrinogen Equivalent Units) µg/L or mg/L |
| Storage temperature | 2–30 °C (36–86 °F) |
| Operating temperature | 15–30 °C (59–86 °F) |
| Measuring range | 190–4000 µg/L FEU (0.190–4 mg/L FEU) Each laboratory should investigate the transferability of the expected values to its own patient population and, if necessary, determine its own reference ranges. |
| Sample type | Whole blood (capillary fingerstick and sodium citrated venous) and sodium citrated plasma samples |
| Sample size | 15 µL |
| Time to result | 6 minutes |

References

1. EMBOL-1 protocol – NCT04737954 www.clinicaltrials.gov/ct2/show/NCT04737954?term=LumiraDx&draw=2&rank=3
2. NICE Guideline: Venous thromboembolism: diagnosis and anticoagulation treatment www.nice.org.uk/guidance/ng158/resources/visual-summary-pdf-8709091453

For more information visit lumiradx.com or contact the LumiraDx Customer Services by email: customerservices@lumiradx.com

This specification sheet and its content are copyright of LumiraDx Group Limited © 2022. All rights reserved, worldwide. Content should be used for the use of the LumiraDx products only and in line with instructions provided. You may not, except with our express written permission, distribute or commercially exploit the content. Nor may you transmit it or store it in any other form of electronic retrieval system other than for the purpose of use of the LumiraDx Instrument or LumiraDx Test Strips. Information provided is subject to change without notice.


Not all products are available in all countries and regions. Please check with your local LumiraDx sales representative or distributor for availability in specific markets.

Not currently available in the USA.

Manufactured by:
LumiraDx UK Ltd,
Dumyat Business Park,
Alloa, FK10 2PB, UK.
Registration Number:
09206123

**Authorized Representative
in the European Union:**
LumiraDx AB
Västra Vägen 5A
16961 Solna,
Sweden

lumiradx.com



S-COM-ART-00328 R3