## SARS-CoV-2 & Flu A/B Quality Controls

For use with the LumiraDx Platform  
For Professional Use Only  
For *In Vitro* Diagnostic Use **IVD**

Product Name	REF
SARS-CoV2 & Flu A/B Quality Controls 2pk	L019080101002
SARS-CoV2 & Flu A/B Quality Controls 2pk	L019080102002

SPEC-34365 R3 ART-01219 R4 Date of Rev 2022-03

The LumiraDx SARS-CoV2 & Flu A/B Quality Controls (hereafter referred to as Quality Controls) are liquid quality controls to be used with the LumiraDx Instrument (hereafter referred to as the Instrument) and the LumiraDx SARS-CoV2 & Flu A/B Test (Penset) referred to as SARS-CoV2 & Flu A/B Test.  
Read these instructions thoroughly before using the Quality Controls.  
Inspect the Quality Controls packaging and contents for damage before use. Report any damage to LumiraDx Customer Services and do not use the kit if any damage is observed to the contents.  
To ensure that you are using the Instrument, the SARS-CoV2 & Flu A/B Test and the Quality Controls correctly, spot the Connect Platform User Manual, SARS-CoV2 & Flu A/B Test Product Insert and the entire product insert. In addition, please watch the LumiraDx Platform training video available at lumiraDx.com.

**Intended Use:**  
The LumiraDx SARS-CoV2 & Flu A/B Quality Controls are intended for use by laboratory professionals/ healthcare professionals for quality control testing of the Instrument. The Quality Controls provide users with assurance that the device is performing with specification.  
The LumiraDx SARS-CoV2 & Flu A/B Quality Control kit should be used whenever local testing guidelines require a QC to be run.  
The Quality Controls are intended for professional use only.

**Summary and Explanation of the Test:**  
The LumiraDx SARS-CoV2 & Flu A/B Quality Controls are liquid quality controls. The controls are specifically formulated and manufactured to ensure performance of the test and are used to verify the user's ability to properly perform the test and interpret the results.  
It is the responsibility of each laboratory or healthcare setting using the LumiraDx SARS-CoV2 & Flu A/B Test to establish an adequate quality assurance program to ensure the performance of the device under its specific locations and conditions of use. Quality control requirements should be followed in conformity with local regulations or accreditation requirements and the user laboratory's standard quality control procedures.

**Warnings and Precautions:**

- For *In vitro* diagnostic use
- This control contains human source material and as with all human based specimens, should be treated as potentially infectious and handled with proper laboratory safety procedures to minimize the risk of transmission of infectious disease.
- All components of this kit can be discarded as biohazard waste according to the local guidelines.
- Refer to the product safety data sheet for risk and safety phrases and disposal information. The product safety data sheet is available for users upon request.
- Requirements of the appropriate licensing or accreditation body should be incorporated into your quality control program.
- Exercise the normal precautions required for handling all laboratory reagents.

**Storage and stability:**

- Store controls between 2°C and 8°C (36- 46°F).
- Unopened controls that are stored between 2°C and 8°C (36 and 46°F) can be used until the expiration date. Do not use kits or components beyond the expiration date given on the label.
- Once opened, the kit has a 30 day expiry.
- Open the Control Vials only when you are performing tests.
- Recap and store the Control Vials in their original container at 2°C and 8°C (36 and 46°F) after use.

**Control Contents:**

- 2 x 0.5ml vial SARS-CoV2 & Flu A/B Positive Control
- 2 x 0.5ml vial SARS-CoV2 & Flu A/B Negative Control
- 40 Transfer pipettes (20µl)
- LumiraDx SARS-CoV2 & Flu A/B Quality Control Pack Insert

**Materials required but not provided with the Control Carton:**

- LumiraDx Instrument
- LumiraDx SARS-CoV2 & Flu A/B Test Strips
- LumiraDx Connect – If connectivity required (refer to LumiraDx Connect User Manual)

**Getting ready to test:**  
You will need the LumiraDx Instrument and the following supplies:  
LumiraDx SARS-CoV2 & Flu A/B Test Strip(s)  
LumiraDx SARS-CoV2 & Flu A/B Positive or Negative Quality Controls  
Transfer pipette

**Pepparing the Quality Controls:**  
The liquid Quality Controls are supplied ready to use.

**Handling the LumiraDx SARS-CoV2 & Flu A/B Test Strips:**  
To ensure that you are using the SARS-CoV2 & Flu A/B Test and the Instrument correctly read the appropriate SARS-CoV2 & Flu A/B Test Strip Product Insert and Platform User Manual.

**UDI** Indicates a carrier that contains unique device identifier information

**CONTROL** Indicates a control material that is intended to verify the performance characteristics of the LumiraDx Instrument. An indication of the level of control is shown in the box, e.g. 1 or 2

**CONTROL-** Indicates a negative control

**CONTROL+** Indicates a positive control

**Bio** Indicates that there are potential biological risks associated with the medical device

**EC REP** Indicates the authorized representative in the European Community/ European Union

**LumiraDx customer service:**  
For product enquiries please contact LumiraDx Customer Services by email: customerservices@lumiraDx.com  
Any adverse results experienced with the use of this product, and/or quality problems should also be reported to LumiraDx Customer Services by email: customerservices@lumiraDx.com or at lumiraDx.com.

**For return policy:**  
If there is a problem with the LumiraDx SARS-CoV2 & Flu A/B Quality Controls Pack you may be asked to return them. Before returning kits, please obtain a return authorization number from LumiraDx Customer Services. This return authorization number must be on the shipping carton for return. For return policy, please contact LumiraDx Customer Services for terms and conditions.

## Procedures/Performing a Test

Consult the LumiraDx Platform User Manual for instructions on how to analyse a Quality Control sample. Open the top pouch of the SARS-CoV2 & Flu A/B Test Strip just before use and insert the Test Strip into the LumiraDx Instrument. Hold the Test Strip by gripping the blue label end with the label facing upward. Do not touch Test Strip Sample Application Area. Do not bend or fold the Test Strip. Do not touch Test Strip contacts. The Instrument will indicate when ready for the sample to be applied.

1. Draw up the Quality Control solution into the Transfer pipette. Squeeze the bulb, drawing up QC material to where the tip of the pipette bases as indicated by the arrow on the diagram. Take care to fill only the tip of the pipette to collect approximately 20µl of QC material.

2. Apply the Quality Control solution to the already inserted SARS-CoV2 & Flu A/B Test Strip. Hold the pipette over the Sample Application Area of the Test Strip and dispense the Quality Control solution. The LumiraDx Instrument will indicate sample is detected with an audible alert (if the Instrument sounds are enabled). The screen of the LumiraDx Instrument will request the user to close the door. Dispose of the pipette.

3. Do not open the door while the test is in progress. The touch-screen will indicate test progress.

4. The Result will appear on the touch-screen within 12 minutes of applying the Quality Control solution and starting the test. The results will be displayed as a PASS or FAIL. QC result for SARS-CoV2 & Flu A/B on the instrument screen.

5. NOTE: If you need to repeat a test, use a new Test Strip.

**Expected results:**  
The Instrument displays the result as Pass or Fail. The result is automatically saved in the memory of the Instrument. The system is working properly, and all handling has been done correctly when the test results obtained are displayed as PASS.

**Symbol glossary:**

Symbol	Meaning
	Temperature limitation
	Manufacturer
<b>IVD</b>	In Vitro Diagnostic Medical Device
<b>REF</b>	Catalogue Number
<b>LOT</b>	Batch code/lot Number
	Use-by Date - indicates the date after which the unopened IVD/Quality Control Material cannot be used
	Refer to Instructions for Use
<b>UK CA</b>	UK conformity assessed under the Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002)
	Indicates the entity importing the medical device into the locale
<b>CE</b>	"CE Mark": This product fulfils the requirements of the European Directive 98/79/EC on in vitro diagnostic medical devices.
<b>UDI</b>	Indicates a carrier that contains unique device identifier information
<b>CONTROL</b>	Indicates a control material that is intended to verify the performance characteristics of the LumiraDx Instrument. An indication of the level of control is shown in the box, e.g. 1 or 2
<b>CONTROL-</b>	Indicates a negative control
<b>CONTROL+</b>	Indicates a positive control
<b>Bio</b>	Indicates that there are potential biological risks associated with the medical device
<b>EC REP</b>	Indicates the authorized representative in the European Community/ European Union

## Warnings - Limited Warranty

LumiraDx SARS-CoV2 & Flu A/B Quality Controls Pack – As per shelf life.  
For the applicable warranty period, LumiraDx warrants that each product shall be (i) of good quality and free of material defects, (ii) function in accordance with the material specifications referenced in the pack insert and (iii) approved by the proper government agencies required for the sale of products for their intended use (the "limited warranty"). If the product fails to meet the requirements of the limited warranty, then on the customer's sole remedy, LumiraDx shall either repair or replace, at LumiraDx's discretion, the LumiraDx SARS-CoV2 & Flu A/B Quality Control Pack. Except for the limited warranty stated in this section, LumiraDx disclaims any and all warranties, express or implied, including but not limited to, any warranty of merchantability, fitness for a particular purpose and non-irritation regarding the product. LumiraDx's maximum liability with any customer claim shall not exceed the net product price paid by the customer. Neither party shall be liable to the other party for special, incidental or consequential damages, including without limitation, loss of business, profits, data or revenue, even if a party receives notice in advance that these kinds of damages might result.  
The Limited Warranty above shall not apply if the customer has subjected the LumiraDx SARS-CoV2 & Flu A/B Test Strips and Control Packs to physical abuse, misuse, abnormal use, use inconsistent with the LumiraDx Platform User Manual, SARS-CoV2 & Flu A/B Test Product Insert or SARS-CoV2 & Flu A/B Quality Control Pack insert head, tampering, unusual physical stress, negligence or accidents. Any warranty claim by Customer pursuant to the Limited Warranty shall be made in writing within the applicable Limited Warranty period.

**Intellectual property:**  
The LumiraDx Instrument, Test Strips and all provided LumiraDx documentation ("Products") are protected by law. The Intellectual Property of the LumiraDx Products remains at LumiraDx. Details of relevant Intellectual Property regarding our products can be found at lumiraDx.com/IP

**Legal notices:**  
Copyright © 2022 LumiraDx UK and affiliates. All rights reserved. LumiraDx and Rame logo are protected trademarks of LumiraDx International LTD. Full details of these and other registrations of LumiraDx can be found at lumiraDx.com/IP. All other trademarks are the property of their respective owners.

**Manufacturer information:**  
LumiraDx UK Ltd, Dumyat Business Park, Alcoa, FK10 2PL UK.  
Registration Number: 09206123

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**CE mark** applies to LumiraDx Instrument, Test Strips, Quality Controls and Connect Hub only.

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Copyright © 2022 LumiraDx UK e affiliazioni. Tutti i diritti riservati. LumiraDx e il logo con la fiamma sono marchi registrati di LumiraDx International LTD. Tutti i dettagli di queste e altre registrazioni di LumiraDx sono disponibili sul sito [lumiraDx.com/IR](http://lumiraDx.com/IR). Tutti gli altri marchi di fabbrica sono di proprietà dei rispettivi titolari.
Copyright © 2022 LumiraDx UK y filiales. Todos los derechos reservados. LumiraDx y el logotipo de la fiamma son marcas comerciales protegidas de LumiraDx International LTD. Para obtener información completa sobre estas y otras registros de LumiraDx, visite [lumiraDx.com/IR](http://lumiraDx.com/IR). Todas las demás marcas comerciales pertenecen a sus respectivos propietarios.