



SARS-CoV-2 Ag Ultra Test Specifications

For *in vitro* Diagnostic Use.

Intended use*

The LumiraDx SARS-CoV-2 Ag Ultra test is an automated rapid microfluidic immunofluorescence assay for use with the LumiraDx Platform, for near-patient testing, intended for the qualitative detection of the nucleocapsid protein antigen from SARS-CoV-2 directly from anterior nasal swab samples collected from individuals suspected of COVID-19 by their healthcare provider within the first twelve days of symptom onset or from individuals without symptoms or other epidemiological reasons to suspect COVID-19.

Test description

The LumiraDx SARS-CoV-2 Ag Ultra test uses SARS-CoV/SARS-CoV-2 specific antibodies in a particle-particle sandwich immunoassay to determine the presence of SARS-CoV-2 Nucleocapsid Protein (NP) antigen present in the test sample.

Clinical performance*

SARS-CoV-2 Ag Ultra test – symptomatic and asymptomatic data

In clinical studies, the LumiraDx SARS-CoV-2 Ag Ultra test demonstrated 97.4%** positive agreement versus RT-PCR in samples to Ct <34 from individuals with symptoms up to and including 12 DSSO.

	Symptomatic	Asymptomatic
PPA	97.4%**	95.7%
NPA	100%	100%

*See SARS-CoV-2 Ag Ultra product insert for full Intended use statement and clinical data

** Ct <34

PPA- POSITIVE PERCENT AGREEMENT; NPA – NEGATIVE PERCENT AGREEMENT; DSSO – DAYS SINCE SYMPTOM ONSET

Built-in quality controls

The LumiraDx Instrument and LumiraDx SARS-CoV-2 Ag Ultra test have several quality control functions integrated to ensure validity of each test run. These checks include:

- Electrical component operation, heater operation, battery charge state, mechanical actuators and sensors and optical system performance.
- Test Strip positioning, optics, and Test Strip expiry.
- Monitoring of Test Strip performance and controls during test runtime.
- The SARS-CoV-2 Ag Ultra test contains an Onboard Quality Control (OBC) assay.

SARS-CoV-2 Ag external quality controls

Positive and Negative Quality Controls are available from LumiraDx to complete Quality Control assessment of the Instrument and SARS-CoV-2 Ag Ultra Test Strips.

Analytical performance

Starting Material Concentration	Estimated LoD	No. Positive/Total	% Positive
1.26 x 10 ⁶ TCID ₅₀ /mL	800 TCID ₅₀ /mL	20/20	100

Cross reactivity

SARS-CoV-2 Ag Ultra test was found not to cross-react with a panel of organisms and viruses including several human coronaviruses. See LumiraDx SARS-CoV-2 Ag Ultra Product Insert for full details.

Specifications

Sample Type	Nasal swabs
Time to Result	5 minutes
Kit Size	24 or 48 Tests
Result Display	Qualitative – Positive or Negative
Storage temperature	2-30 °C (36-86 °F)
Operating temperature	15-30 °C (59-86 °F)
Interferences	See LumiraDx SARS-CoV-2 Ag Ultra Product Insert for details
Onboard Control	Onboard Quality Control (OBC) assay and sample processing control
Quality Control Material	Positive and Negative external liquid controls
Nasal Collection Swabs	Sterile nasal collection swabs available with certain product codes

Swabs

Individually packaged sterile nasal collection swabs are included in certain kits and should be used where provided.

For information on swabs that have been validated for use with the LumiraDx SARS-CoV-2 Ag Ultra test please visit lumiradx.com and consult the technical bulletin.

For more information visit lumiradx.com or contact the LumiraDx Customer Services by email: customerservices@lumiradx.com or Phone: 0080058647239


Copyright © 2022 LumiraDx UK LTD. All rights reserved worldwide.

LumiraDx and Flame logo are trademarks of LumiraDx International LTD. Full details of these and other registrations of LumiraDx can be found at lumiradx.com/IP. All other trademarks are the property of their respective owners.

Content should be used for the use of the LumiraDx products only and in line with instructions provided. You may not, except with our express written permission, distribute or commercially exploit the content. Nor may you transmit it or store it in any other form of electronic retrieval system other than for the purpose of use of the LumiraDx Instrument or LumiraDx Test Strips. Information provided is subject to change without notice.

Product is not available in all countries and regions. Please check with your local LumiraDx sales representative or distributor for availability in specific markets.

lumiradx.com



Manufactured by:
LumiraDx UK Ltd,
Dumyat Business Park,
Alloa,
FK10 2PB, UK.

Registration Number:
09206123

**Authorized Representative
in the European Union:**
LumiraDx AB
Västra Vägen 5A
16961 Solna, Sweden

S-COM-ART-02502 R1