 LumiraDx™ SARS-CoV-2 Ag Pool

For Professional Use Only  
For *In Vitro* Diagnostic Use Only

SPEC-33561 R3 ART-00897 R6

<b>ENGLISH</b>	01 - 29
<b>FRANÇAIS</b>	31 - 59
<b>DEUTSCH</b>	61 - 89
<b>ITALIANO</b>	91 - 119
<b>ESPAÑOL</b>	121 - 149
<b>NEDERLANDS</b>	151 - 178

# LumiraDx™ SARS-CoV-2 Ag Pool

For Professional Use Only  
For *In Vitro* Diagnostic Use Only

Product name	REF	
LumiraDx SARS-CoV-2 Ag Pool Test	L016000201048	48

SPEC-33561 R3    ART-00897 R6

## LumiraDx SARS-CoV-2 Ag Pool

The LumiraDx Severe Acute Respiratory Syndrome (SARS) CoV-2 Antigen (Ag) Pool Test Strips (hereafter referred to as Test Strips) are to be used with the LumiraDx Platform. The LumiraDx Platform is a point of care system for professional use which is used for *in vitro* diagnostic tests. It comprises a portable LumiraDx Instrument and a LumiraDx Test Strip for the required test. This test is for **HEALTHCARE PROFESSIONAL USE ONLY** and allows users to perform tests using small sample volumes and to view results quickly on the Instrument touchscreen.

### Intended use:

The LumiraDx SARS-CoV-2 Ag Pool Test is a rapid microfluidic immunofluorescence assay for use with the LumiraDx Platform intended for the qualitative detection of the nucleocapsid protein antigen to SARS-CoV-2 in 1 to 5 individual samples from professionally supervised & self-collected nasal swab samples or professionally collected nasal or nasopharyngeal swab samples which are then pooled for testing. Samples should be collected from 1 to 5 individuals suspected of COVID-19 infection within the first twelve days of symptom onset or from asymptomatic individuals. The Test aids in the diagnosis of current SARS-CoV-2 infection by detection of SARS-CoV-2 antigen.

Positive results indicate the presence of viral antigens from infective virus in one of the pooled samples, but clinical correlation with each individual's history and other diagnostic information is necessary to confirm infection status. If the pool is positive, then each individual should be re-tested as a separate and newly collected individual sample before reporting a result.

Negative results do not rule out SARS-CoV-2 infection and should be considered in the context of an individual's recent exposures, history and presence of clinical signs and symptoms consistent with COVID-19.

Results should not be used as the sole basis for treatment or case management decisions, including infection control decisions.

The LumiraDx SARS-CoV-2 Ag Pool Test is intended for use by individuals trained in point of care settings and proficient in performing tests using the LumiraDx Platform.

**Caution:** For *in vitro* diagnostic use.

Before you start testing, if you are new to the LumiraDx Instrument and LumiraDx Platform, you must read the LumiraDx Platform User Manual and this entire Product Insert. In addition, please watch the LumiraDx Platform Training Video available at [kc.lumiradx.com](http://kc.lumiradx.com).

## Summary and explanation of the Test:

The World Health Organization (WHO) have named the disease caused by SARS-CoV-2 virus as coronavirus 2019 or COVID-19<sup>1</sup>. The most common symptoms of COVID-19 are fever, tiredness, and dry cough. Some patients may have aches and pains, nasal congestion, headache, conjunctivitis, sore throat, diarrhoea, loss of taste or smell, or a rash on skin or discoloration of fingers or toes. These symptoms are usually mild and begin gradually. Some people become infected but do not develop any symptoms and do not feel unwell. However, the disease can develop rapidly and have high morbidity in certain populations, especially those with underlying health conditions. The disease can spread from person to person through small droplets from the nose or mouth which are spread when a person with COVID-19 coughs or exhales. Most estimates of the incubation period for COVID-19 range from 2-14 days<sup>2</sup>.

The use of a LumiraDx SARS-CoV-2 Ag Pool Test will enable the physician to verify infection quickly among a pool of 1 to 5 patients, begin proper treatment and initiate isolation precautions helping prevent further spread of infection.

## Principle of the assay:

The LumiraDx SAR-CoV-2 Ag Pool Test is a fluorescence immunoassay device designed to detect the presence of the nucleocapsid protein antigen from SARS-CoV-2 in nasal swab and nasopharyngeal swab samples. The Test Strip is single use.

The test procedure involves collecting 1 to 5 nasal swab or nasopharyngeal swab samples using a recommended swab which are then eluted into a vial containing Extraction Buffer. A single drop of the sample in Extraction Buffer is added to the Test Strip using the vial dropper cap provided. The LumiraDx Instrument is programmed to perform the test protocol using the dried reagents contained within the Test Strip. The test result is determined from the amount of fluorescence the Instrument detects within the measurement zone of the Test Strip. The concentration of the analyte in the sample is proportional to the fluorescence detected. The results are displayed on the Instrument touchscreen within 12 minutes from the addition of the sample.

In sample pooling, the collected nasal or nasopharyngeal swab samples are identified and up to five individual samples are then combined into one pool by eluting each individual sample sequentially into a single vial containing Extraction Buffer. Attach the dropper cap, invert the sample 5x, and squeeze to place a single drop of the pooled sample from the Extraction Buffer to the Test Strip. The LumiraDx Instrument is programmed to perform the Pooled Test Protocol within 12 minutes from the addition of the pooled sample. If the pool is positive, then each individual should be re-tested as a separate and newly collected individual sample. If the pool is negative, then each constituent sample is reported as negative.

## Materials provided:

- LumiraDx SARS-CoV-2 Ag Pool Test
  - LumiraDx Test Strips packed individually in sealed desiccant foil pouches.
  - Extraction Buffer Vial
  - Dropper Lid
  - LumiraDx Test Product Insert
  - RFID (Radio frequency ID) Tag held inside the Test Strip carton

## Materials required but not provided:

- Swabs and collection tubes
  - Please refer to list of approved swabs on website
  - Use a dry, clean, and sterile tube
- LumiraDx Instrument
- LumiraDx SARS-CoV-2 Ag Pool Test Quick Reference Instructions
- LumiraDx SARS-CoV-2 Ag Quality Controls (as required to meet local and organisational compliance)

### **Materials not required but available:**

- LumiraDx Connect if connectivity required (refer to LumiraDx Connect User Manual)
- LumiraDx Barcode Scanner
- LumiraDx How to Use Platform video
- LumiraDx How to Use an Antigen Test video

"How to" video's are accessible at [www.lumiradx.com](http://www.lumiradx.com)

### **Warnings and precautions:**

- For *in vitro* diagnostic use only
- Do not open the Test Strip until ready for immediate use.
- Discard and do not use any damaged or dropped Test Strips or other materials.
- Inadequate or inappropriate sample collection, storage, and transport can result in incorrect results.
- All samples in a single pool must be of the same sample type (nasal or nasopharyngeal)
- The Test cannot be visually interpreted; the LumiraDx Instrument must be used to generate results.
- Do not use the kit components beyond the expiration date
- Do not reuse any kit components.
- Samples must be processed as indicated in the Sample Extraction and Performing a Test sections of this Product Insert. Failure to follow the instructions for use can result in inaccurate results.
- All components of this kit should be discarded as biohazard waste according to local regulations and procedures.
- Refer to the product safety data sheet for risk and safety phrases and disposal information. The product safety data sheet is available via our website at [lumiradx.com](http://lumiradx.com).
- Exercise the normal precautions required for handling all laboratory reagents. Wear protective clothing such as laboratory coats, disposable gloves, and eye protection when samples are collected and evaluated.
- Proper laboratory safety techniques should be followed at all times when working with SARS-CoV-2 patient samples. Patient swabs, used Test Strips, and used Extraction Buffer vials may be potentially infectious. Proper handling and disposal methods should be established by the laboratory in accordance with local regulations and procedures.
- For additional information on safety, handling, and disposal of the components within this kit, please refer to the Safety Data Sheet (SDS) located at [lumiradx.com](http://lumiradx.com).

### **Storing the Test Kit:**

Store the Test Strips in their original carton. You can store the Test Strips at a temperature between 2°C and 30°C (36°F and 86°F). Avoid freezing or storing in any area that could exceed 30°C. When stored properly, the Test Strips can be used until the expiration date printed on the Test Strip foil pouch and the Test Strip carton. Discard the Test Strips if they are passed the expiration date.

### **Handling the Test Strips:**

When you are ready to perform a test, open the Test Strip carton, take out a Test Strip, and remove it from the foil pouch. Hold the Test Strip by gripping the blue label end with the label facing upward. Do not touch Test Strip Sample Application Area. Do not bend or fold the Test Strip. Do not touch Test Strip contacts. After removing the Test Strip from the foil pouch, it should be used immediately. Do not use the Test Strip if there are any visible signs of damage to the foil pouch such as tears or holes.

## Sample material:

The following samples can be used with the LumiraDx SARS-CoV-2 Ag Pool Test Strip:

- Nasal Swab Samples (NS)
- Nasopharyngeal Swab Samples (NP)

## The Test device contains:

- Rabbit and Mouse Monoclonal antibodies
- Fluorescent Latex particles
- Magnetic particles
- Buffer and Stabilising agents

## Preparing the Instrument to perform a Test:

Power on the Instrument by pressing the power button at the rear of the Instrument. You will hear the Instrument powering on, and the display will be a blank black screen for several seconds before starting up. If the screen is just dimmed tap the touch-screen to wake up the Instrument. Check that "Pooled Test" is available on the Instrument home screen. If it is not, enable "Pooled Test" in the Instrument settings menu.

If performing the SARS-CoV-2 Ag Pool Test with 2-5 patients, select "Pooled Test" from the Platform Home Screen. If performing the SARS-CoV-2 Ag Pool Test with 1 patient, select "Patient Test" from the Platform Home Screen. Please do not insert the Test Strip until prompted to do so.

Refer to the section on **Performing a Test** in this Product Insert for information on how to test a sample. The LumiraDx SARS-CoV-2 Ag Pool Test Quick Reference Instructions (QRI) provide an illustrated step-by-step procedure on how to run a Test. Operate the LumiraDx Platform at room temperature between 15°C and 30°C (59°F and 86°F) and 10% - 75% relative humidity.

The Instrument will prompt to install the Lot Calibration File when inserting a new Test Strip Lot for the first time. Once installed, the Instrument will have all the information required to process the test, and any future tests from the same Lot of Test Strips.

## Lot Calibration File installation

Lot Calibration Files are required to provide the Instrument with the information needed to perform diagnostic tests. This only needs to be completed once for each Test Strip Lot. The Instrument will prompt to install the Lot Calibration File when inserting a new Test Strip Lot. To install, locate the ((•)) symbol on the instrument and touch it to the same symbol on the back of the Test Strip Carton to install.

### RFID strip code reader

Locate ((•)) symbol on Instrument.

### Installation

Touch back of Test Strip Carton ((•)) symbol to install.

The Instrument will sound and a confirmation message will be displayed.

When indicated by the touchscreen, open the foil pouch just before use and insert the LumiraDx Test Strip into the LumiraDx Instrument. The Instrument will indicate when it is ready for the sample to be applied. The LumiraDx SARS-CoV-2 Ag Pool Test results should be evaluated by a Healthcare Professional in the context of all available clinical and laboratory data.


### Instructions for sample collection:

When collecting any type of sample, follow universal collection precautions and guidelines according to your organization. For collection of nasal swabs and nasopharyngeal swabs, follow appropriate Swab Collection Guidelines and swab manufacturers' recommendations. Users should be trained in appropriate sample collection and handling procedures.


The steps that follow apply to a nasal swab and a nasopharyngeal swab. For information on recommended swabs to use with the LumiraDx SARS-CoV-2 Ag Pool Test, please see the "SARS-CoV-2 Ag Technical Bulletin – Swabs" available at Lumiradx.com.

### Sampling from a nasal swab:


Individual samples should be collected and placed into a dry, clean, and sterile tube separately. The steps that follow apply to a nasal swab.

1. Tilt patient's head back 70°


2. A swab sample is needed from both nostrils and this is taken using the same swab. Remove sterile swab from the swab packet. Hold the swab by the shaft, while gently rotating the swab, insert swab less than one inch into the nostril until resistance is met at Turbinates. (Turbinates are the small structures inside the nose).


3. Rotate the swab several times against the nasal wall. Remove and repeat this process by using the same swab into the second nostril. Place swab in dry, clean, and sterile tube or process the swab directly in the extraction buffer vial as per instructions for sample extraction of pooled samples outlined below.

### Sampling from a nasopharyngeal swab:


Individual samples should be collected and placed into a dry, clean, and sterile tube separately. The steps that follow apply to a nasopharyngeal swab.

1. Tilt patient's head back 70°

2. Remove sterile swab from the packet. Hold the swab firmly by the shaft between the fingers and insert Swab into nostril. The swab should reach a depth equal to distance from nostrils to outer opening of the ear. Leave Swab in place for several seconds to absorb secretions.

3. Slowly remove the swab while rotating it. Remove and then place the swab in a dry, clean, and sterile tube or process the swab directly in the extraction buffer vial as per instructions for sample extraction of pooled samples outlined below.

**NOTE:** After patient swabbing, process the Swab(s) in the Extraction Buffer Vial as soon as possible. Freshly collected patient samples should be processed in the Extraction Buffer vial within 1 hour of collection. If the sample is not processed immediately, the patient swab may be stored at room temperature (15°C - 30°C/59°F - 86°F) for up to 1 hour prior to testing in a dry, clean and sterile tube.

Do not place the swab back into the swab packaging sleeve after sample collection.


### Sample pooling for SARS-CoV-2 testing

Users may consider not using pooled testing or changing the pool size to test fewer than 5 pooled samples based on a variety of factors, including SARS-CoV-2 pre-test prevalence. Pools of 1 to 5 swab samples may be tested using the SARS-CoV-2 Ag Pool Test.


Process pooled samples as described in the instructions below for sample extraction of pooled samples.

#### Instructions for sample extraction of pooled samples:


1 to 5 individual swabs can be eluted sequentially into a single Extraction Buffer vial.

1. **Remove the seal or blue screw cap** from the top of the Extraction Buffer Vial containing the Extraction Buffer.


2. **Place and soak the Patient Swab** in the Extraction Buffer for 10 seconds and then stir well by rotating the swab against the side of the vial 5 times.


3. **Squeeze and remove** the Patient Swab while squeezing the middle of the Extraction Buffer Vial to remove the liquid from the swab. Discard the swab in biohazard waste.

**NOTE: For testing of a pooled samples of 2-5 samples, repeat steps 2 and 3 sequentially** for up to 4 more swabs. Up to and including a total of 5 samples can be placed into the same Extraction Buffer Vial.

4. **Firmly attach the clear or purple Dropper Lid** to the top of the Extraction Buffer Vial. The extracted sample must be used within 5 hours of preparation when stored at room temperature. Extracted nasal swab or nasopharyngeal swab samples may be frozen at  $-80^{\circ}\text{C}$  and used up to 5 days after freezing.


5. Gently **invert the Extraction Buffer Vial** 5 times just before applying the sample to the Test Strip.

## Performing a Test (refer to the Quick Reference Instructions to make sure that your Instrument has been prepared before starting this step)


1. **Squeeze and apply the extracted sample from the Extraction Buffer Vial** onto the Sample Application Area of the inserted Test Strip. To do this gently press the sides of the extraction buffer vial until **one whole drop** is visible and allow it to touch the Sample Application Area of the Test Strip. The sample will then be drawn by capillary action into the Test Strip. When the sample is detected the Instrument will sound (if sounds are enabled) and a confirmation message will be displayed. The touchscreen of the LumiraDx Instrument will request the user to **immediately close the door (Note: you have 10 seconds only to close the door)**.
2. **Do not add more than one drop of sample.** Do not open the door while the test is in progress. The touchscreen will indicate test progress.
3. **The result** will appear on the Instrument touchscreen within 12 minutes of applying the sample and starting the test. For pooled sample testing with 2-5 samples, the results will be displayed as a **"Pooled Positive"** or **"Pooled Negative" result SARS-CoV-2 Ag** on the Instrument screen. (see Fig 1 and Fig 2). For pooled sample testing with 1 sample, the results will be displayed as a **positive or negative result SARS-CoV-2 Ag** on the Instrument screen.
4. **If you need to retest**, you will use a new Test Strip. Use the same Extraction Buffer Vial and repeat the test. The extracted sample must be used within 5 hours of preparation when stored at room temperature. Extracted nasal swab or nasopharyngeal swab samples may be frozen at -80°C and used up to 5 days after freezing.
5. **Dispose** of the swab, Extraction Buffer Vial and Test Strip in the appropriate clinical waste.
6. **Disinfection** of the Instrument with LumiraDx approved materials is recommended if contamination is suspected. Details of approved disinfecting materials is available at lumiradx.com. Allow the Instrument to air dry before testing the next sample. The disinfectant should remain in contact for at least 1 minute.

### Result interpretation for individual test:

The results will be displayed on the Instrument screen – **examples are shown below of result screen display for pools of 2-5:**

**Fig 1: Negative result for SARS-CoV-2 Ag Pool Test of 2-5 patients**

**Fig 2: Positive result for SARS-CoV-2 Ag Pool Test of 2-5 patients**

**NOTE: A negative result, from patients with symptoms onset beyond twelve days, should be treated as presumptive and confirmation with a molecular assay, if necessary, for patient management, may be performed.**

**Negative**— Negative results from pooled sample testing require no further testing of individuals within the pool and each constituent sample is reported as negative. If the individual's clinical signs and symptoms are inconsistent with a negative result and if results are necessary for individual sample management, then the individual should be considered for individual testing.


**Positive**— A positive pool means that one or more of the individuals tested in that pool may be positive for SARS-CoV-2 antigen. Individuals included within a positive sample pool with the LumiraDx SARS-CoV-2 Ag Pool Test should (1) be recalled for testing with individual sample collection or (2) seek follow up care with their physician or healthcare provider for additional testing.

## Error Messages

If an issue occurs, a message will be displayed on the Instrument touch-screen. Alert messages include useful information and are highlighted by an orange banner. Error messages also include a symbol. All messages will contain a description of the Instrument status or error and an instruction. Error messages contain an identifying code that may be used for further troubleshooting purposes. Refer to the LumiraDx Platform User Manual if an error message is displayed on the LumiraDx Instrument touch-screen and contact LumiraDx Customer Services on [customerservices@lumiradx.com](mailto:customerservices@lumiradx.com).

### Example of an error screen:

If the On Board Control (OBC) fails, an error message will be shown and no test result will be returned. Follow the on screen instructions to dispose of the Test Strip and start a new test. If the problem persists, contact Customer Services.

## Built-in controls:

The instrument reads the 2D bar code on each Test Strip and can identify if the strip has exceeded the expiry date for use, and if the strip Lot Calibration file has not yet been loaded, at which point it will request it.

The LumiraDx Instrument and LumiraDx SARS-CoV-2 Ag Pool Test Strips have several quality control functions integrated to ensure validity of each test run. These checks ensure that the volume of sample added is sufficient and the assay sequence of the Test Strip is as expected. The checks also ensure that the Test Strip has not been damaged or used previously. If these checks are not verified, the test run will be rejected, and an error message displayed on the Instrument touchscreen.

The LumiraDx Instrument ensures the quality of test results obtained through the following features:

- Automated checks of the correct functioning of the Instrument at power on and during operation.
- This includes electrical component operation, heater operation, battery charge state, mechanical actuators and sensors and optical system performance.
- Monitoring of Test Strip performance and controls during test run time.
- Ability to perform Quality Control Tests using LumiraDx Quality Control solutions to meet regulatory compliance requirements.

## External Quality Controls:

External liquid SARS-CoV-2 Ag Quality Controls are available from LumiraDx and may be used to demonstrate that the Test is functioning properly by demonstrating the expected Quality Control results and correct test performance by the operator.

External Quality Control requirements should be established in accordance with local, state, and federal regulations or accreditations requirements. It is recommended that external control testing be performed with each new operator and before using a new lot or shipment of the LumiraDx SARS-CoV-2 Ag Pool Test. Refer to the LumiraDx SARS-CoV-2 Ag Quality Controls pack insert available at [www.lumiradx.com](http://www.lumiradx.com) for detailed instructions.

LumiraDx SARS-CoV-2 Ag Quality Controls are purchased separately.

If the LumiraDx SARS-CoV-2 Ag Quality Controls do not perform as expected, repeat the QC Test and if the problems persist, do not report patient results and contact LumiraDx Customer Services.

## Cleaning and disinfection:

Cleaning and disinfection of the Instrument should follow and be performed according to established site protocols and schedules.

To clean the Instrument wipe the external surfaces with a soft, slightly damp cloth when it appears visibly dirty.

It is recommended to clean and disinfect the Instrument at least once per day if in use, or if contamination is suspected. Details of LumiraDx approved disinfectant materials can be found at [lumiradx.com](http://lumiradx.com). Allow the Instrument to air dry before testing the next sample. The disinfectant should remain in contact for at least 1 minute.

Excessive liquid may damage the Instrument. It is important for the protection of the Instrument that exposure to excess moisture is prevented. All disinfection cloths and/or wipes should only be slightly damp, with any excess liquid being manually removed from the cloth before use.

**Avoid USB ports and power inlet. Do not spray or pour solution directly onto the Instrument. Do not put any objects or cleaning materials into the Test Strip slot.**

## Limitations:

- This test detects both viable (live) and non-viable, SARS-CoV and SARS-CoV-2. Test performance depends on the amount of virus (antigen) in the sample and may or may not correlate with viral culture results performed on the same sample.
- Failure to follow the instructions for use may adversely affect test performance and/or invalidate the test result.
- Test results should be considered in the context of all available clinical and diagnostic information, including patient history and other test results.
- Positive test results do not differentiate between SARS-CoV and SARS-CoV-2.
- Negative test results are not intended to rule in other non-SARS viral or bacterial infections.
- Negative results, from patients with symptom onset beyond twelve days, should be treated as presumptive and confirmation with a molecular assay, if necessary for patient management, may be performed.
- If the differentiation of specific SARS viruses and strains is needed, additional testing, in consultation with state or local public health departments, is required.
- When considering pooling strategies, the appropriateness of the pooling strategy should be considered based on the positivity rate in the testing population, efficiency of the pooling workflow, and pooled test specificity.
- Users should test samples as quickly as possible after sample collection.

- Samples extracted into the same Extraction Buffer should all be of the same sample type (nasal or nasopharyngeal).
- Extracted nasal samples or nasopharyngeal samples may be frozen at -80°C and used up to 5 days after freezing.
- Swab samples and Extraction Buffer must be at room temperature before testing.
- Positive test results do not rule out co-infection with other pathogens
- A false negative result may occur if the level of viral antigen in a sample is below the detection limit of the test or if the sample was collected inappropriately, therefore a negative test result does not rule out the possibility of SARS-CoV-2 infection.
- The amount of antigen in a sample may decrease as the duration of illness increases. Samples collected after 12 days are more likely to be negative compared to RT-PCR.
- The contents of this kit are for qualitative detection of SARS-CoV-2 antigens from nasal swab and nasopharyngeal swab samples only.
- For information on swabs have been validated for use with the LumiraDx SARS-CoV-2 Ag Pool Test please visit [lumiradx.com](http://lumiradx.com).
- Only dry, clean, and sterile tubes have been validated for use with the LumiraDx SARS-CoV-2 Ag Pool Test.
- Frozen swabs have not been validated for individual or pooled testing.

## **Clinical performance**

### **SARS-CoV-2 Ag Test with 1 Nasal Sample**

The performance of the LumiraDx SARS-CoV-2 Ag Test was established with 257 direct nasal swabs prospectively collected from individual subjects during the 2020 COVID-19 pandemic. Samples were collected from sequentially enrolled subjects who presented with symptoms of COVID-19 (159) or key workers (98) at increased risk of infection. No positive results were observed from patients without symptoms or beyond 12 days of symptom onset. Dual nasal swabs were simultaneously collected and then randomly allocated to testing with the LumiraDx test or the Roche Cobas 6800. Samples were collected from 6 sites across the United States (5) and United Kingdom (1), including four sites in which minimally trained operators collected and tested fresh samples.

Swabs were collected and extracted into the LumiraDx Extraction Buffer without transport media. Samples were tested fresh or frozen within 1h of collection and stored until tested. Samples were thawed and sequentially tested according to the Product Insert, with operators blinded to the PCR result. The performance of the LumiraDx SARS-CoV-2 Ag Pool Test was compared to the results from nasal swabs collected into 3mL universal transport medium (UTM) and tested with the Roche Cobas 6800 PCR method.

## Patient demographics

Patient demographics (gender, age, time elapsed since onset of symptoms) are available for the 257 samples used in the study. The table below shows the positive results broken down by age of the patient:

Age	LumiraDx SARS-CoV-2 Ag (n = 81)		
	Total #	Positive	Prevalence
≤ 5 years	13	0	0.0%
6 to 21 years	29	6	20.7%
22 to 59 years	200	70	35.0%
≥ 60 years	15	5	33.3%

Positive results broken down by days since symptom onset:

Days since symptom onset	Cumulative RT-PCR Positive(+)	Cumulative LumiraDx Positive(+)	PPA	95% Confidence interval	
0	6	6	100.0%	61.0%	100.0%
1	12	12	100.0%	75.8%	100.0%
2	28	28	100.0%	87.9%	100.0%
3	37	37	100.0%	90.6%	100.0%
4	55	54	98.2%	90.4%	99.7%
5	61	60	98.4%	91.3%	99.7%
6	67	66	98.5%	92.0%	99.7%
7	73	72	98.6%	92.6%	99.8%
8	75	74	98.7%	92.8%	99.8%
9	75	74	98.7%	92.8%	99.8%
10	77	76	98.7%	93.0%	99.8%
11	80	79	98.8%	93.3%	99.8%
12	83	81	97.6%	91.6%	99.3%

Final data analysis is presented below:

Reference RT-PCR Assay							95% Wilson Score CI		
							LCI	UCI	
LumiraDx SARS-CoV-2 Ag Test		POS	NEG	Total	PPA	97.6%	91.6%	99.3%	
	POS	81	6	87	NPA	96.6%	92.7%	98.4%	
	NEG	2	168	170	PPV	93.1%	85.8%	96.8%	
	TOTAL	83	174	257	NPV	98.8%	95.8%	99.7%	
						Prevalence	32.3%	26.9%	38.2%
						OPA	96.9%	94.0%	98.4%

100% positive agreement was achieved with samples with cycle threshold values (Ct) < 33 (\*Roche Cobas 6800 SARS-CoV-2 reference method).

Evidence<sup>3</sup> suggests that patients with Ct values above 33-34 are no longer contagious.

Published literature shows that RT-PCR Ct values from symptomatic and asymptomatic subjects are similar<sup>4,5</sup>.

The LumiraDx Antigen Test detects the nucleocapsid protein antigen in extracted nasal samples in both symptomatic and asymptomatic subjects.

**PPA** - Positive Percent Agreement (Sensitivity)

**NPA** - Negative Percent Agreement (Specificity)

**PPV** - Positive Predictive Value

**NPV** - Negative Predictive Value

**OPA** - Overall Percent Agreement

**CI** - Confidence Interval

**LCI** - Lower Confidence Interval

**UCI** - Upper Confidence Interval

### SARS-CoV-2 Ag Test with 1 Nasopharyngeal Sample

The performance of the SARS-CoV-2 Ag Test was established with 255 nasopharyngeal swabs prospectively collected from individual subjects during the 2020 COVID pandemic. Subjects were presenting with symptoms of COVID-19 being screened for infection. Samples were collected from 6 sites across the United States. Swabs were collected and extracted into the LumiraDx Extraction Buffer. Samples were tested fresh within 1h of collection and tested according to the Product Insert. The performance of the LumiraDx SARS-COV-2 Ag Test was compared to the results from nasopharyngeal samples collected into 3ml universal transport medium (UTM) and tested with an EUA authorized PCR method.

## Patient Demographics

Patient demographics (gender, age, time elapsed since onset of symptoms) are available for the 255 samples used in the study. The table below shows the positive results broken down by age of the patient:

Age	LumiraDx SARS-CoV-2 Ag (n = 39)		
	Total #	Positive	Prevalence
≤ 5 years	22	0	0.0%
6 to 21 years	59	9	15.3%
22 to 59 years	150	28	18.7%
≥ 60 years	24	2	8.3%

Positive and negative results broken down by days since symptom onset:

Days since symptom onset	Cumulative RT-PCR Positive(+)	LumiraDx Positive(+)	PPA	95% Confidence Interval		NPA	95% Confidence Interval	
				LCI	UCI		LCI	UCI
0	2	2	100.0%	34.2%	100.0%	100.0%	75.8%	100.0%
1	6	6	100.0%	61.0%	100.0%	100.0%	93.4%	100.0%
2	9	9	100.0%	70.1%	100.0%	100.0%	96.2%	100.0%
3	17	17	100.0%	81.6%	100.0%	98.6%	94.9%	99.6%
4	22	22	100.0%	85.1%	100.0%	98.8%	95.7%	99.7%
5	23	23	100.0%	85.7%	100.0%	98.4%	95.3%	99.4%
6	26	26	100.0%	87.1%	100.0%	98.5%	95.6%	99.5%
7	34	34	100.0%	89.8%	100.0%	98.5%	95.7%	99.5%
8	36	36	100.0%	90.4%	100.0%	98.6%	95.8%	99.5%
9	36	36	100.0%	90.4%	100.0%	98.6%	95.9%	99.5%
10	39	38	97.4%	86.8%	99.5%	98.1%	95.2%	99.3%
11	40	39	97.5%	87.1%	99.6%	97.7%	94.6%	99.0%
12	40	39	97.5%	87.1%	99.6%	97.7%	94.7%	99.0%

Final Data analysis is presented below:

Reference RT-PCR Assay							95% Wilson Score CI	
							LCI	UCI
LumiraDx SARS-CoV-2 Ag Test		POS	NEG	Total	PPA	97.5%	87.1%	99.6%
	POS	39	5	44	NPA	97.7%	94.7%	99.0%
	NEG	1	210	211	PPV	88.6%	76.0%	95.0%
	TOTAL	40	215	255	NPV	99.5%	97.4%	99.9%
					Prevalence	15.7%	11.7%	20.7%
					OPA	97.6%	95.0%	98.9%

100% positive agreement was achieved with samples with cycle threshold values (ct) < 33 (\*Roche Cobas 6800 SARS-CoV-2 reference method).

Evidence<sup>3</sup> suggests that patients with Ct values above 33-34 are no longer contagious. Published literature shows that RT-PCR Ct values from symptomatic and asymptomatic subjects are similar<sup>4,5</sup>.

### Clinical Performance – Nasal and Nasopharyngeal Swabs Symptomatic and Asymptomatic Subjects

Clinical studies were continued to include asymptomatic and symptomatic positive individuals across sites in the United Kingdom, the United States and Norway\*. Studies were conducted with individual subjects.

The performance of the LumiraDx SARS-CoV-2 Ag Test was established with direct nasal swabs prospectively collected from asymptomatic and symptomatic individual subjects during the 2020 COVID-19 pandemic. Samples were collected from 782 subjects in the US, UK and Norway (176 positive by PCR). Dual nasal swabs were simultaneously collected and then randomly allocated to testing with the LumiraDx test or the RT-PCR reference. Swabs were collected and extracted into the LumiraDx Extraction Buffer and tested fresh or frozen within 1h of collection and stored until tested. The performance of the LumiraDx SARS-CoV-2 Ag Test was compared to the results from nasal swabs or nasopharyngeal swabs (NP) collected into 3ml universal transport medium (UTM) and tested with RT-PCR. Ct range = 15-38.

The performance of the LumiraDx SARS-CoV-2 Ag Test was established with nasopharyngeal swabs prospectively collected from asymptomatic and symptomatic individual subjects during the 2020 COVID-19 pandemic. Samples were collected from 676 subjects in the US and Norway (109 positive by PCR). Dual NP swabs were collected and then randomly allocated to testing with the LumiraDx Test or the RT-PCR reference. Swabs were collected and extracted into the LumiraDx Extraction Buffer and used fresh within 1h of collection. The performance of the LumiraDx SARS-CoV-2 Ag Test was compared to the results from NP swabs collected into 3ml universal transport medium (UTM) and tested with RT-PCR.

\* SKUP (Scandinavian evaluation of laboratory equipment for point of care testing) independent evaluation of the LumiraDx SARS-CoV-2 Ag test. [www.skup.org](http://www.skup.org).

Grouping	Nasal Swabs					NP Swabs				
	N	PPA	95% CI	NPA	95% CI	N	PPA	95% CI	NPA	95% CI
<b>All Subjects</b>	<b>782</b>	<b>93.8%</b>	<b>89.2-96.5</b>	<b>98.8%</b>	<b>97.6-99.4</b>	<b>676</b>	<b>92.7%</b>	<b>86.2-96.2</b>	<b>98.1%</b>	<b>96.6-98.9</b>
DSSO* ≤ 12	473	95.0%	90.4-97.4	98.4%	96.3-99.3	484	93.9%	87.4-97.2	97.9%	96.0-98.9
DSSO ≤ 10	465	95.4%	90.9-97.8	98.4%	96.3-99.3	479	93.9%	87.3-97.2	98.2%	96.3-99.1
DSSO ≤ 7	451	95.3%	90.6-97.7	98.3%	96.2-99.3	463	94.7%	88.1-97.7	98.4%	96.5-99.3
Ct < 33 (all)**	146	97.9%	94.1-99.3	-	-	103	97.1%	91.8-99.0	-	-
Ct < 30 (all)	129	99.2%	95.7-99.9	-	-	93	97.8%	92.5-99.4	-	-
Ct < 25 (all)	91	100%	95.9-100	-	-	73	98.6%	92.6-99.8	-	-
Asymptomatic	309	82.4%	59.0-93.8	99.3%	97.5-99.8	192	80.0%	49.0-94.3	98.4%	95.3-99.4
Ct < 33 (asymptomatic)	13	100%	77.2-100	-	-	7	100%	64.6-100	-	-

\* DSSO = Days Since Symptom Onset

\*\* Ct < (X) (all) = includes all subjects in the study (symptomatic and asymptomatic)

### Clinical Performance – Asymptomatic Subjects

The performance of the SARS-CoV-2 Ag Test was established with 222 anterior nasal swabs prospectively collected from individual asymptomatic subjects between June 2020 and March 2021. Samples were collected from 5 sites across the United States. Swabs were collected and extracted into the LumiraDx Extraction Buffer. Samples were tested fresh according to the Product Insert. The performance of the LumiraDx SARS-COV-2 Ag Test was compared to the results from paired anterior nasal swab samples collected into 3ml universal transport medium (UTM) and tested with an EUA authorized PCR method. Data analysis is presented below:

### Patient Demographics

Patient demographics (age, prevalence) are available for the 222 samples used in the study. The table below shows the positive results broken down by age of the patient:

Age	LumiraDx SARS-CoV-2 Ag		
	Total #	Positive	Prevalence
≤ 5 years	0	0	0.0%
6 to 21 years	49	7	14.3%
22 to 59 years	145	11	7.6%
≥ 60 years	28	5	17.9%

Reference RT-PCR Assay							95% Wilson Score CI	
							LCI	UCI
LumiraDx SARS-CoV-2 Ag Test		POS	NEG	Total	PPA	82.1%	64.4%	92.1%
	POS	23	0	23	NPA	100.0%	98.1%	100.0%
	NEG	5	194	199	PPV	100%	85.7%	100%
	TOTAL	28	194	222	NPV	97.5%	94.3%	98.9%
					Prevalence	12.6%	8.9%	17.6%
					OPA % Agreement	97.7%	94.8%	99.0%

PPA - Positive Percent Agreement (Sensitivity)

NPA – Negative Percent Agreement (Specificity)

PPV – Positive Predictive Value

NPV - Negative Predictive Value

OPA – Overall Percent Agreement

At Ct < 33, the PPA was 95.8% and 100% PPA for Ct < 30

#### SARS-CoV-2 Ag Pool Test with 5 Nasal Samples

The performance of the SARS-CoV-2 Ag Pool Test was established during the 2020 COVID pandemic with 30 positive pools prospectively collected from 150 individual subjects, and 29 negative pools prospectively collected from 145 samples. Subjects were either presenting with symptoms of COVID-19 or key workers being screened for infection. Samples (anterior nares) were collected by a healthcare provider or by supervised on site self-collection from 1 site each in the United States and United Kingdom. Swabs were collected and extracted into the LumiraDx extraction buffer. Samples were extracted and tested fresh within 1h of collection. One positive swab was included in a randomized sequence with 4 negative swabs for the positive pools. For the negative pools, 5 negative swabs were extracted. The performance of the LumiraDx SARS-CoV-2 Ag Pool Test was compared to the results from nasal swabs collected into the LumiraDx extraction buffer and tested with the LumiraDx SARS-COV-2 Ag Test.

Final data analysis is presented below.

Classification of all evaluable pools:

		Single Swab		TOTAL
		POS	NEG	
Pooled Test	POS	30	1	31
	NEG	0	28	28
TOTAL		30	29	59

Performance measures, and 95% confidence intervals, as calculated with the Wilson Score method.

Measure	Estimate	95% Confidence Interval	
PPA	100.0%	88.6%	100.0%
NPA	96.6%	82.8%	99.4%
PPV	96.8%	83.8%	99.4%
NPV	100.0%	87.9%	100.0%
Prevalence	50.8%	38.4%	63.2%
OPA (% Agreement)	98.3%	91.0%	99.7%

PPA - Positive Percent Agreement (Sensitivity)

NPA - Negative Percent Agreement (Specificity)

PPV - Positive Predictive Value

NPV- Negative Predictive Value

OPA - Overall Percent Agreement

#### Analytical performance:

##### Limit of Detection - LoD (Analytical sensitivity):

Limit of Detection (LoD) studies determined the lowest detectable concentration of SARS-CoV-2 at which 100% of all (true positive) replicates test positive. The LoD for the LumiraDx SARS-CoV-2 Ag Test was established using limiting dilutions of gamma-irradiated SARS-CoV-2 (BEI Resources NR-52287). The NR-52287 is a preparation of SARS-Related Coronavirus 2 (SARS-CoV2), isolate USA WA1/2020, that has been inactivated by gamma-irradiation at  $5 \times 10^6$  RADs. The material was supplied frozen at a concentration of  $2.8 \times 10^5$  TCID<sub>50</sub>/mL.

##### **Limit of Detection (LoD) screening**

An initial LoD screening study was performed using a 5-fold serial dilutions (six dilutions in total) of the gamma-irradiated virus made in pooled negative human nasal matrix starting at a test concentration of  $2 \times 10^4$  TCID<sub>50</sub>/mL (as shown in table below) and processed for each study as described above. These dilutions were tested in triplicate. The lowest concentration at which all (3 out of 3 replicates) were positive was chosen for LoD Range finding. This was 32 TCID<sub>50</sub>/mL.

SARS-CoV-2 tested (TCID <sub>50</sub> /mL)	Test result
20000	3/3 positive
4000	3/3 positive
800	3/3 positive
160	3/3 positive
32	3/3 positive
6.4	0/3 positive

### Limit of Detection range finding

Using the 32 TCID<sub>50</sub>/mL concentration, the LoD was further refined using a 2-fold dilution series (four dilutions in total) of the gamma-irradiated SARS-CoV-2 virus made in pooled negative human nasal matrix. These dilutions were tested in triplicate. The lowest concentration at which all (3 out of 3 replicates) were positive was treated as the tentative LoD for the LumiraDx SARS-CoV-2 Ag Test. This was 32 TCID<sub>50</sub>/mL.

SARS-CoV-2 tested (TCID <sub>50</sub> /mL)	Test result
32	3/3 positive
16	0/3 positive
8	1/3 positive
4	0/3 positive

### Limit of Detection (LoD) confirmation

The LoD of the LumiraDx SARS-CoV-2 Ag Test was then confirmed by testing 20 replicates with concentrations at the tentative Limit of Detection. The final LoD of the LumiraDx SARS-CoV-2 Ag Test was determined to be the lowest concentration resulting in positive detection of twenty (20) out of twenty (20) replicates.

Starting Material Concentration	Estimated LoD	No. Positive/Total	% Positive
2.8 x 10 <sup>5</sup> TCID <sub>50</sub> /mL	32 TCID <sub>50</sub> /mL	20/20	100

The Limit of Detection studies were repeated for the SARS-CoV-2 Pool Test with one positive swab and 4 negative swabs and demonstrated comparable results to the SARS-CoV-2 Ag Test when run across 3 test strip batches.

### Cross-reactivity (analytical specificity) and microbial interference studies with 1 nasal sample

Cross-reactivity and interference of the LumiraDx SARS-CoV-2 Ag Test was evaluated by testing a panel of related pathogens, high prevalence disease agents and normal or pathogenic flora including various microorganisms and viruses and negative matrix that are reasonably likely to be encountered in the clinical sample and could potentially cross-react or interfere with the LumiraDx SARS CoV-2 Ag Test. Each organism and virus were tested in the absence or presence of heat inactivated SARS-CoV-2 at 3 x LoD and results are shown below.

Microorganism	Source	Concentration	Cross-reactivity (Yes/No)	Interference (Yes/No)
Human coronavirus 229E	Zeptomatrix	1 x 10 <sup>5</sup> PFU/mL	No (3/3 negative)	No (3/3 positive)
Human coronavirus OC43	Zeptomatrix	1 x 10 <sup>5</sup> PFU/mL	No (3/3 negative)	No (19/20 positive)
Human coronavirus NL63	Zeptomatrix	9.87 x 10 <sup>3</sup> PFU/mL	No (3/3 negative)	No (3/3 positive)
MERS coronavirus	Zeptomatrix	7930 PFU/mL	No (2/2 negative)	No (3/3 positive)

Microorganism	Source	Concentration	Cross-reactivity (Yes/No)	Interference (Yes/No)
Adenovirus (e.g. C1 Ad. 71)	Zeptomatrix	1 x 10 <sup>5</sup> PFU/mL	No (3/3 negative)	No (3/3 positive)
Human Metapneumovirus (hMPV)	Zeptomatrix	1 x 10 <sup>5</sup> PFU/mL	No (3/3 negative)	No (3/3 positive)
Parainfluenza virus Type 1	Zeptomatrix	1 x 10 <sup>5</sup> PFU/mL	No (3/3 negative)	No (3/3 positive)
Parainfluenza virus Type 2	Zeptomatrix	1 x 10 <sup>5</sup> PFU/mL	No (3/3 negative)	No (3/3 positive)
Parainfluenza virus Type 3	Zeptomatrix	1 x 10 <sup>5</sup> PFU/mL	No (3/3 negative)	No (3/3 positive)
Parainfluenza virus Type 4a	Zeptomatrix	1 x 10 <sup>5</sup> PFU/mL	No (3/3 negative)	No (3/3 positive)
Influenza A H3N2 (Wisconsin/67/05)	Zeptomatrix	8.82 x 10 <sup>4</sup> PFU/mL	No (3/3 negative)	No (3/3 positive)
Influenza A H1N1	Zeptomatrix	1 x 10 <sup>5</sup> PFU/mL	No (3/3 negative)	No (3/3 positive)
Influenza B (Malaysia/2506/04)	Zeptomatrix	2.92 x 10 <sup>4</sup> PFU/mL	No (3/3 negative)	No (19/20 positive)
Enterovirus	Zeptomatrix	1 x 10 <sup>5</sup> PFU/mL	No (3/3 negative)	No (3/3 positive)
Respiratory syncytial virus	Zeptomatrix	1 x 10 <sup>5</sup> PFU/mL	No (3/3 negative)	No (3/3 positive)
Rhinovirus	Zeptomatrix	4.17 x 10 <sup>5</sup> PFU/mL	No (3/3 negative)	No (3/3 positive)
<i>Haemophilus influenzae</i>	Zeptomatrix	1 x 10 <sup>6</sup> CFU/mL	No (3/3 negative)	No (3/3 positive)
<i>Streptococcus pneumoniae</i>	Zeptomatrix	1 x 10 <sup>6</sup> CFU/mL	No (3/3 negative)	No (3/3 positive)
<i>Streptococcus pyogenes</i>	Zeptomatrix	1 x 10 <sup>6</sup> CFU/mL	No (3/3 negative)	No (3/3 positive)
<i>Candida albicans</i>	Zeptomatrix	1 x 10 <sup>6</sup> CFU/mL	No (3/3 negative)	No (3/3 positive)
Pooled human nasal wash	LumiraDx	14% v/v	No (3/3 negative)	No (3/3 positive)
<i>Bordetella pertussis</i>	Zeptomatrix	1 x 10 <sup>6</sup> CFU/mL	No (3/3 negative)	No (3/3 positive)
<i>Mycoplasma pneumoniae</i>	ATCC	1 x 10 <sup>6</sup> CFU/mL	No (3/3 negative)	No (3/3 positive)
<i>Chlamydia pneumoniae</i>	ATCC	1 x 10 <sup>6</sup> CFU/mL	No (3/3 negative)	No (3/3 positive)
<i>Legionella pneumophila</i>	Zeptomatrix	1 x 10 <sup>6</sup> CFU/mL	No (3/3 negative)	No (3/3 positive)
<i>Mycobacterium tuberculosis</i>	Zeptomatrix	1 x 10 <sup>6</sup> CFU/mL	No (3/3 negative)	No (3/3 positive)

Microorganism	Source	Concentration	Cross-reactivity (Yes/No)	Interference (Yes/No)
<i>Pneumocystis jirovecii</i>	Zeptomatrix	1 x 10 <sup>6</sup> CFU/mL	No (3/3 negative)	No (3/3 positive)
<i>Pseudomonas Aeruginosa</i>	Zeptomatrix	1 x 10 <sup>6</sup> CFU/mL	No (3/3 negative)	No (3/3 positive)
<i>Staphylococcus Epidermidis</i>	Zeptomatrix	1 x 10 <sup>6</sup> CFU/mL	No (3/3 negative)	No (3/3 positive)
<i>Streptococcus Salivarius</i>	Zeptomatrix	1 x 10 <sup>6</sup> CFU/mL	No (3/3 negative)	No (3/3 positive)
<i>Staphylococcus aureus</i>	ATCC	1 x 10 <sup>6</sup> CFU/mL	No (3/3 negative)	No (3/3 positive)

To estimate the likelihood of cross-reactivity with SARS-CoV-2 of organisms that were not available for wet testing, *in silico* analysis using the Basic Local Alignment Search Tool (BLAST) managed by the National Center for Biotechnology Information (NCBI) was used to assess the degree of protein sequence homology.

- For Human Coronavirus HKU1, homology exists between the SARS-CoV-2 nucleocapsid protein and Human Coronavirus HKU1. BLAST results showed 30 sequence IDs, all nucleocapsid protein, showing homology. Sequence ID AGW27840.1 had the highest alignment score and was found to be 39.1% homologous across 76% of the sequences, this is relatively low but cross-reactivity cannot be fully ruled out.
- For SARS-Coronavirus, high homology exists between the SARS-CoV-2 nucleocapsid protein and SARS-Coronavirus. BLAST results showed 68 sequence IDs, mostly nucleocapsid protein, showing homology. Sequence ID AAR87518.1, had the highest alignment score isolated from a human patient and was found to be 90.76% homologous across 100% of the sequence. This is high and cross-reactivity is likely.
- For MERS-Coronavirus, high homology exists between the SARS-CoV-2 nucleocapsid protein and MERS-Coronavirus. BLAST results showed at least 114 sequence IDs, mostly nucleocapsid protein, showing homology. Sequence IDs AHY61344.1 and AWH65950.1, had the highest alignment scores isolated from a human patient and were found to be 49.4% and 50.3% homologous across 88% of the sequence. Whilst this potentially represents moderate cross-reactivity testing of the MERS virus at 7930 PFU/mL showed no reactivity (see table above).

#### **Cross-reactivity (analytical specificity) and microbial interference studies with 5 nasal samples**

Higher risk cross-reactant and microbial interferents were additionally tested using the LumiraDx SARS-CoV-2 Ag Pool Test. The following were tested:

- Human coronavirus OC43
- Adenovirus (e.g. C1 Ad. 71)
- Parainfluenza virus 1
- Parainfluenza virus 2
- Parainfluenza virus 3
- Parainfluenza virus 4b
- Influenza B
- *Pneumocystis jirovecii* (PJP)

Cross-reactivity and interference of the LumiraDx SARS-CoV-2 Ag Pool Test was evaluated by testing a panel of high-risk pathogens that are reasonably likely to be encountered in the clinical sample and could potentially cross-react or interfere with the LumiraDx SARS CoV-2 Ag Pool Test. Each organism or virus testing concentration was spiked to each of the 5 pooled swabs and tested in the absence or presence of heat inactivated SARS-CoV-2 at 3 x LoD. The results are shown below.

Microorganism	Source	Concentration	Cross-reactivity (Yes/No)	Interference (Yes/No)
Human coronavirus OC43	Zeptomatrix	1 x 10 <sup>5</sup> PFU/mL	No (3/3 negative)	No (3/3 positive)
Adenovirus (e.g. C1 Ad. 71)	Zeptomatrix	1 x 10 <sup>5</sup> PFU/mL	No (3/3 negative)	No (3/3 positive)
Parainfluenza virus Type 1	Zeptomatrix	1 x 10 <sup>5</sup> PFU/mL	No (3/3 negative)	No (3/3 positive)
Parainfluenza virus Type 2	Zeptomatrix	1 x 10 <sup>5</sup> PFU/mL	No (3/3 negative)	No (3/3 positive)
Parainfluenza virus Type 3	Zeptomatrix	1 x 10 <sup>5</sup> PFU/mL	No (3/3 negative)	No (3/3 positive)
Parainfluenza virus Type 4b	Zeptomatrix	1 x 10 <sup>5</sup> PFU/mL	No (3/3 negative)	No (3/3 positive)
Influenza B	Zeptomatrix	2.72 x 10 <sup>4</sup> PFU/mL	No (3/3 negative)	No (3/3 positive)
<i>Pneumocystis jirovecii</i>	Zeptomatrix	1 x 10 <sup>6</sup> CFU/mL	No (3/3 negative)	No (3/3 positive)

#### Endogenous and Exogenous interference studies with 1 nasal sample

A study was performed to demonstrate that potentially interfering substances that may be found in the upper respiratory tract in symptomatic subjects (including over the counter medications) do not cross-react or interfere with the detection of SARS-CoV-2 in the LumiraDx SARS-CoV-2 Ag Test. Each substance was tested in triplicate in the absence or presence of SARS-CoV-2 at 3 x LoD. The results are shown below.

Interfering substance	Concentration	Interference (Yes/No)
Benzocaine**	150 mg/dL	No (3/3 Negative, 3/3 Positive)
Blood (human)*	5%	No (3/3 Negative, 3/3 Positive)
Mucin*	5 mg/mL	No (3/3 Negative, 3/3 Positive)
Naso GEL (NeilMed)**	5% v/v	No (3/3 Negative, 3/3 Positive)
CVS Nasal Drops (phenylephrine)**	15% v/v	No (3/3 Negative, 3/3 Positive)

<b>Interfering substance</b>	<b>Concentration</b>	<b>Interference (Yes/No)</b>
Afrin (Oxymetazoline)**	15% v/v	No (3/3 Negative, 3/3 Positive)
CVS Nasal Spray (Cromolyn)**	15% v/v	No (3/3 Negative, 3/3 Positive)
Zicam Cold Remedy**	5% v/v	No (3/3 Negative, 3/3 Positive)
Homeopathic (Alkaloi)**	10% v/v	No (3/3 Negative, 3/3 Positive)
Sore Throat Phenol Spray**	15% v/v	No (3/3 Negative, 3/3 Positive)
Tobramycin**	3.3 mg/dL	No (3/3 Negative, 3/3 Positive)
Mupirocin**	0.15 mg/dL	No (3/3 Negative, 3/3 Positive)
Fluticasone**	0.000126 mg/dL	No (5/5 Negative, 4/4 Positive)
Tamiflu (Osetamivir phosphate)**	500 mg/dL	No (3/3 Negative, 3/3 Positive)
Budesonide**	0.00063 mg/dL	No (3/3 Negative, 3/3 Positive)
Biotin**	0.35 mg/dL	No (3/3 Negative, 3/3 Positive)
Methanol**	150 mg/dL	No (19/20 Negative, 3/3 Positive)
Acetylsalicylic Acid**	3 mg/dL	No (3/3 Negative, 3/3 Positive)
Diphenhydramine**	0.0774 mg/dL	No (3/3 Negative, 3/3 Positive)
Dextromethorphan**	0.00156 mg/dL	No (19/20 Negative, 3/3 Positive)
Dexamethasone**	1.2 mg/dL	No (3/3 Negative, 3/3 Positive)
Mucinex**	5%	No (3/3 Negative, 3/3 Positive)

\* Endogenous substances

\*\* Exogenous substances

## Endogenous and Exogenous interference studies with 5 nasal samples

Higher risk interferents were additionally tested using the LumiraDx SARS-CoV-2 Ag Pool Test. The following were tested:

- Acetylsalicylic acid
- Dexamethasone

A study was performed to demonstrate that potentially interfering substances that may be found in the upper respiratory tract in symptomatic subjects do not cross-react or interfere with the detection of SARS-CoV-2 in the LumiraDx SARS-CoV-2 Ag Pool Test. Each testing concentration was spiked to each of the 5 pooled swabs and was tested in triplicate in the absence or presence of SARS-CoV-2 at 3 x LoD. The results are shown below.

Interfering substance	Concentration	Interference (Yes/No)
Acetylsalicylic Acid**	3 mg/dL	No (6/6 Negative, 6/6 Positive)
Dexamethasone**	1.2 mg/dL	No (6/6 Negative, 6/6 Positive)

\* Endogenous substances

\*\* Exogenous substances

## High dose hook effect

High Dose Hook Effect studies determine the level at which false negative results can be seen when very high levels of target are present in a tested sample. To determine if the LumiraDx SARS-CoV-2 Ag Test suffers from any high dose hook effect, increasing concentrations of gamma-irradiated SARS-CoV-2 virus (BEI Resources NR-52287) were tested up to a concentration of  $1.4 \times 10^5$  TCID<sub>50</sub>/mL. In this study, the starting material was spiked into a volume of pooled human nasal matrix obtained from healthy donors and confirmed negative for SARS-CoV-2. At each dilution, 50  $\mu$ L samples were added to swabs and the swabs processed for testing on the LumiraDx SARS-CoV-2 Ag Test as per the Product Insert using the procedure appropriate for patient nasal swab samples.

No impact on test performance or high dose hook effect was observed up to  $1.4 \times 10^5$  TCID<sub>50</sub>/mL of gamma-irradiated SARS-CoV-2 with the LumiraDx SARS-CoV-2 Ag Test.

Test Dilution	Concentration (TCID <sub>50</sub> /mL)	Mean Signal (ADC Units)
1	0	495
2	62.5	26100.6
3	250	63013.8
4	1000	83451.8
5	$1.4 \times 10^5$	86220



## Variants of Concern

LumiraDx actively monitors new mutations in the SARS-CoV-2 viral genome as they arise. The reactivity of the LumiraDx SARS CoV-2 Ag test will be assessed against all variants of concern as they arise. The up to date results of this testing program can be found on our SARS-CoV-2 Variants Technical Bulletin on our website [lumiradx.com](http://lumiradx.com).

## References:

1. World Health Organization [www.who.int](http://www.who.int)
2. Centers for Disease Control and Prevention [www.cdc.gov](http://www.cdc.gov)
3. La Scola B., Le Bideau M., Andreani J., Hoang V.T., Grimaldier C., Colson P. Viral RNA load as determined by cell culture as a management tool for discharge of SARS-CoV-2 patients from infectious disease wards. *Eur J Clin Microbiol Infect Dis*. 2020:1–3.
4. Kissler, SM., et al., Viral dynamics of SARS-CoV-2 infection and the predictive value of repeat testing <https://www.medrxiv.org/content/10.1101/2020.10.21.20217042v1>
5. Lavenzo, E., et al., Suppression of a SARS-CoV-2 outbreak in the Italian municipality of Vo'. *Nature* 584, 2020: 425–429

## Symbols glossary:

Symbol	Meaning
	Temperature limitation
	Manufacturer
	<i>In vitro</i> diagnostic medical device
	Catalogue Number
	Lot Number
	Use-by Date - indicates the date after which the unopened IVD/Quality Control Material cannot be used
	Refer to instructions for use
	Do not re-use
	Indicates the total number of IVD tests that can be performed with the IVD medical device.
	Indicates a carrier that contains Unique Device Identifier information
	Indicates the authorized representative in the European Community/European Union
	Indicates the presence of the Radio Frequency Identification (RFID) reader/tag.
	For near patient testing
	UK conformity assessed under the Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002)
	Importer
	"CE Mark". This product fulfils the requirements of the European Directive 98/79/EC on <i>in vitro</i> diagnostic medical devices.

### **LumiraDx customer services:**

For product enquiries please contact LumiraDx Customer Services at [customerservices@lumiradx.com](mailto:customerservices@lumiradx.com) or find telephone contact details at [lumiradx.com](http://lumiradx.com).

Any adverse results experienced with the use of this product, and/or quality problems should also be reported to LumiraDx Customer Services by email: [customerservices@lumiradx.com](mailto:customerservices@lumiradx.com) or at [lumiradx.com](http://lumiradx.com).

### **For return policy:**

If there is a problem with the **LumiraDx SARS-CoV-2 Ag Pool Test Strips** you may be asked to return them. Before returning tests please obtain a return authorization number from LumiraDx Customer Services. This return authorization number must be on the shipping carton for return. For ordinary returns following purchase, please contact LumiraDx Customer Services for terms and conditions.

### **Limited warranty:**

**LumiraDx SARS-CoV-2 Ag Pool Test Strips** – As per shelf life.

Unused strips must be stored according to the required storage conditions as printed in this product insert and they can be used only up to the expiry date printed on the Test Strip pouch and Test Strip box. For the applicable warranty period, LumiraDx warrants that each product shall be (i) of good quality and free of material defects, (ii) function in accordance with the material specifications referenced in the product insert, and (iii) approved by the proper governmental agencies required for the sale of products for their intended use (the "limited warranty"). If the product fails to meet the requirements of the limited warranty, then as customer's sole remedy, LumiraDx shall either repair or replace, at LumiraDx's discretion, the Test Strips. Except for the limited warranty stated in this section, LumiraDx disclaims any and all warranties, express or implied, including but not limited to, any warranty of merchantability, fitness for a particular purpose and non-infringement regarding the product. LumiraDx's maximum liability with any customer claim shall not exceed the net product price paid by the customer. Neither party shall be liable to the other party for special, incidental or consequential damages, including, without limitation, loss of business, profits, data or revenue, even if a party receives notice in advance that these kinds of damages might result. The Limited Warranty above shall not apply if the customer has subjected the LumiraDx SARS-CoV-2 Ag Pool Test to physical abuse, misuse, abnormal use, use inconsistent with the LumiraDx Platform User Manual or Product Insert, fraud, tampering, unusual physical stress, negligence or accidents. Any warranty claim by Customer pursuant to the Limited Warranty shall be made in writing within the applicable Limited Warranty period.

### **Intellectual property:**


The LumiraDx Instrument, Test Strips and all provided LumiraDx documentation ("Products") are protected by law. The Intellectual Property of the LumiraDx Products remains at LumiraDx. Details of relevant Intellectual Property regarding our products can be found at [lumiradx.com/IP](http://lumiradx.com/IP).

### **Legal notices:**

Copyright © 2022 LumiraDx UK and affiliates. All rights reserved. LumiraDx and Flame logo are protected trademarks of LumiraDx International LTD. Full details of these and other registrations of LumiraDx can be found at [lumiradx.com/IP](http://lumiradx.com/IP). All other trademarks are the property of their respective owners.

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CE mark applies to LumiraDx Instrument, Test Strips, Quality Controls and Connect Hub only.

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
### **Literaturangaben:**

1. Weltgesundheitsorganisation [www.who.int](http://www.who.int)
2. Centers for Disease Control and Prevention [www.cdc.gov](http://www.cdc.gov)
3. La Scola B., Le Bideau M., Andreani J., Hoang V.T., Grimaldier C., Colson P. Viral RNA load as determined by cell culture as a management tool for discharge of SARS-CoV-2 patients from infectious disease wards. *Eur J Clin Microbiol Infect Dis.* 2020:1–3.
4. Kissler, SM., et al., Viral dynamics of SARS-CoV-2 infection and the predictive value of repeat testing <https://www.medrxiv.org/content/10.1101/2020.10.21.20217042v1>
5. Lavenzo, E., et al., Suppression of a SARS-CoV-2 outbreak in the Italian municipality of Vo'. *Nature* 584, 2020: 425–429

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
1. Organizzazione Mondiale della Sanità [www.who.int](http://www.who.int)
2. Centers for Disease Control and Prevention [www.cdc.gov](http://www.cdc.gov)
3. La Scola B., Le Bideau M., Andreani J., Hoang V.T., Grimaldier C., Colson P. Viral RNA load as determined by cell culture as a management tool for discharge of SARS-CoV-2 patients from infectious disease wards. *Eur J Clin Microbiol Infect Dis.* 2020:1–3.
4. Kissler, SM., et al., Viral dynamics of SARS-CoV-2 infection and the predictive value of repeat testing <https://www.medrxiv.org/content/10.1101/2020.10.21.20217042v1>
5. Lavenzo, E., et al., Suppression of a SARS-CoV-2 outbreak in the Italian municipality of Vo'. *Nature* 584, 2020: 425–429

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## Literatuur:

1. World Health Organization [www.who.int](http://www.who.int)
2. Centers for Disease Control and Prevention [www.cdc.gov](http://www.cdc.gov)
3. La Scola B., Le Bideau M., Andreani J., Hoang V.T., Grimaldier C., Colson P. Viral RNA load as determined by cell culture as a management tool for discharge of SARS-CoV-2 patients from infectious disease wards. *Eur J Clin Microbiol Infect Dis.* 2020:1–3.
4. Kissler, SM., et al., Viral dynamics of SARS-CoV-2 infection and the predictive value of repeat testing <https://www.medrxiv.org/content/10.1101/2020.10.21.20217042v1>
5. Lavenzo, E., et al., Suppression of a SARS-CoV-2 outbreak in the Italian municipality of Vo'. *Nature* 584, 2020: 425–429

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Copyright © 2022 LumiraDx UK und Tochtergesellschaften. Alle Rechte vorbehalten. LumiraDx und das Flammen-Logo sind geschützte Marken von LumiraDx International LTD. Vollständige Angaben zu diesen und anderen Eintragungen von LumiraDx sind unter [lumiradx.com/IP](http://lumiradx.com/IP) zu finden. Alle anderen Marken sind Eigentum des jeweiligen Inhabers.
Copyright © 2022 LumiraDx UK e affiliate. Tutti i diritti riservati. LumiraDx e il logo con la fiamma sono marchi protetti di LumiraDx International LTD. Tutti i dettagli di queste e altre registrazioni di LumiraDx sono disponibili sul sito [lumiradx.com/IP](http://lumiradx.com/IP). Tutti gli altri marchi di fabbrica sono di proprietà dei rispettivi titolari.
Copyright © 2022 LumiraDx UK y filiales. Todos los derechos reservados. LumiraDx y el logotipo de la llama son marcas comerciales protegidas de LumiraDx International LTD. Para obtener información completa sobre estos y otros registros de LumiraDx, visite [lumiradx.com/IP](http://lumiradx.com/IP). Todas las demás marcas comerciales pertenecen a sus respectivos propietarios.
Copyright © 2022 LumiraDx UK en gelieerde ondernemingen. Alle rechten voorbehouden. LumiraDx en het vlamlogo zijn beschermde handelsmerken van LumiraDx International LTD. Volledige gegevens over deze en andere gedeponeerde handelsmerken van LumiraDx vindt u op [lumiradx.com/IP](http://lumiradx.com/IP). Alle overige handelsmerken zijn eigendom van de desbetreffende merkhouders.
LumiraDx UK Ltd, Dumyat Business Park, Alloa, FK10 2PB, VK.
LumiraDx AB, Västra Vågen 5A, 16961 Solna, Zweden